

ABX Minotrol Retic



LOT

RX 091

CONTROL



(Exp.) 2011-11-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		120 RETIC	DX120		120 RETIC	DX120		120 RETIC	DX120				
Méthode Manuelle / Manual RETIC	%	1.3	1.3		± 0.7	5.6	5.6		± 2.2	11.1	11.1		± 3.3
RETIC	%	1.39	1.46		± 0.73	5.26	5.31		± 1.33	10.70	10.80		± 2.20
RETIC	10 ⁹ /µl	0.048	0.050		± 0.025	0.180	0.180		± 0.045	0.365	0.366		± 0.073
MFI MRI *		31.7	32.4		± 8.0	43.3	43.9		± 10.0	45.3	46.3		± 10.0
GR RBC	10 ⁹ /µl	3.45	3.43		± 0.30	3.43	3.39		± 0.30	3.41	3.39		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C

Ref: FORM-0382 Rev.07

* For research use only in Canada and USA