

ABX Minotrol Retic



LOT RX 052
Rev 1

CONTROL

(Exp.) **2012-07-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX															
		CONTROL			1	TOLERANCES TOLERANCE	CONTROL			2	TOLERANCES TOLERANCE	CONTROL			3	TOLERANCES TOLERANCE	
		PENTRA	PENTRA	PENTRA			PENTRA	PENTRA	PENTRA			PENTRA	PENTRA	PENTRA			PENTRA
		120 RETIC	DX120	DX Nexus			120 RETIC	DX120	DX Nexus			120 RETIC	DX120	DX Nexus			
Méthode Manuelle / Manual																	
RETIC	%	1.5	1.5	1.5	± 0.8	5.6	5.6	5.6	± 2.2	11.1	11.1	11.1	± 3.3				
RETIC	%	1.79	1.83	1.83	± 0.92	5.09	5.06	5.06	± 1.27	10.80	10.90	10.90	± 2.20				
RETIC	10 ⁶ /µl	0.061	0.063	0.063	± 0.032	0.173	0.175	0.175	± 0.044	0.366	0.373	0.373	± 0.075				
MFI MRI *		27.2	28.6	28.6	± 8.0	39.5	40.0	40.0	± 10.0	43.4	44.3	44.3	± 10.0				
GR RBC	10 ⁶ /µl	3.41	3.47	3.47	± 0.30	3.40	3.45	3.45	± 0.30	3.39	3.42	3.42	± 0.30				
PIC RANGE		25.5	25.5	25.5	± 4.5	25.5	25.5	25.5	± 4.5	25.5	25.5	25.5	± 4.5				

9930059-A
Ref: TEMP-0382 Rev.08

* For research use only in Canada and USA