

ABX Minotrol Retic



LOT RX 112
Rev 1

CONTROL

(Exp.) **2013-01-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX															
		CONTROL			1	TOLERANCES TOLERANCE	CONTROL			2	TOLERANCES TOLERANCE	CONTROL			3	TOLERANCES TOLERANCE	
		PENTRA	PENTRA	PENTRA			PENTRA	PENTRA	PENTRA			PENTRA	PENTRA	PENTRA			PENTRA
		120 RETIC	DX120	DX Nexus			120 RETIC	DX120	DX Nexus			120 RETIC	DX120	DX Nexus			
Méthode Manuelle / Manual																	
RETIC	%	1.3	1.3	1.3	± 0.7	5.2	5.2	5.2	± 2.1	10.7	10.7	10.7	± 3.2				
RETIC	%	1.51	1.74	1.74	± 0.87	4.51	4.86	4.86	± 1.22	10.20	10.50	10.50	± 2.10				
RETIC	10 ⁶ /µl	0.052	0.060	0.060	± 0.030	0.154	0.165	0.165	± 0.041	0.344	0.350	0.350	± 0.070				
MFI MRI *		26.7	24.5	24.5	± 8.0	37.9	38.8	38.8	± 10.0	42.0	44.0	44.0	± 10.0				
GR RBC	10 ⁶ /µl	3.48	3.45	3.45	± 0.30	3.41	3.40	3.40	± 0.30	3.37	3.33	3.33	± 0.30				
PIC RANGE		25.5	25.5	25.5	± 4.5	25.5	25.5	25.5	± 4.5	25.5	25.5	25.5	± 4.5				

9930059-A
Ref: TEMP-0382 Rev. 10

* For research use only in Canada and USA