

ABX Minotrol Retic



LOT RX 423
Rev 1

CONTROL

(Exp.) 2020-07-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | | | |
| | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | |
| Méthode Manuelle / Manual RETIC | % | 1.8 | 1.8 | 1.8 | ± 1.0 | 5.4 | 5.4 | 5.4 | ± 2.2 | 10.6 | 10.6 | 10.6 | ± 3.5 |
| RETIC | % | 1.52 | 1.63 | 1.63 | ± 0.82 | 5.01 | 5.11 | 5.11 | ± 1.28 | 10.60 | 10.70 | 10.70 | ± 2.20 |
| RETIC | 10 ⁶ /µl | 0.053 | 0.058 | 0.058 | ± 0.029 | 0.175 | 0.182 | 0.182 | ± 0.046 | 0.373 | 0.383 | 0.383 | ± 0.079 |
| MFI MRI * | % | 30.0 | 29.4 | 29.4 | ± 8.0 | 45.2 | 45.2 | 45.2 | ± 10.0 | 48.7 | 49.9 | 49.9 | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.48 | 3.55 | 3.55 | ± 0.30 | 3.50 | 3.57 | 3.57 | ± 0.30 | 3.52 | 3.58 | 3.58 | ± 0.30 |
| PIC RANGE | - | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 |

Ref: TEMP-0382 Rev.37 FRONT / RECTO 1219930083-A

* For research use only in Canada and USA

ABX Minotrol Retic



LOT RX 423
Rev 1

CONTROL

(Exp.) 2020-07-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | | |
| Méthode Manuelle / Manual RETIC | % | 1.8 | | | ± 1.0 | 5.4 | | | ± 2.2 | 10.6 | | | ± 3.5 |
| RETIC | % | 1.75 | | | ± 0.88 | 5.02 | | | ± 1.26 | 10.50 | | | ± 2.10 |
| RETIC | 10 ⁶ /µl | 0.059 | | | ± 0.030 | 0.171 | | | ± 0.043 | 0.358 | | | ± 0.072 |
| MFI MRI * | % | 26.8 | | | ± 8.0 | 41.4 | | | ± 10.0 | 45.9 | | | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.38 | | | ± 0.30 | 3.40 | | | ± 0.30 | 3.41 | | | ± 0.30 |
| PIC RANGE | - | 8 | | | ± 4 | 8 | | | ± 4 | 8 | | | ± 4 |

Ref: TEMP-0382 Rev.37 BACK / VERSO 1219930083-A

* For research use only in Canada and USA