

ABX Minotrol Retic



LOT RX 435
Rev 1

CONTROL

(Exp.) 2022-07-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|-----------|----------|---------|-------------------------|-----------|---------|----------|-------------------------|-----------|--|--|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | | | | |
| | | DX120 | DX Nexus | DX120 | | DX Nexus | DX120 | DX Nexus | | | | | |
| Méthode Manuelle / Manual RETIC | % | 1.3 | 1.3 | ± 0.7 | 5.2 | 5.2 | ± 2.1 | 10.1 | 10.1 | ± 3.0 | | | |
| RETIC | % | 1.59 | 1.59 | ± 0.80 | 4.97 | 4.97 | ± 1.24 | 10.20 | 10.20 | ± 2.04 | | | |
| RETIC | 10 ⁶ /µl | 0.055 | 0.055 | ± 0.028 | 0.170 | 0.170 | ± 0.043 | 0.346 | 0.346 | ± 0.069 | | | |
| MFI MRI * | % | 28.9 | 28.9 | ± 8.0 | 47.6 | 47.6 | ± 10.0 | 51.6 | 51.6 | ± 12.0 | | | |
| GR RBC | 10 ⁶ /µl | 3.47 | 3.47 | ± 0.30 | 3.43 | 3.43 | ± 0.30 | 3.39 | 3.39 | ± 0.30 | | | |
| PIC RANGE | - | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | ± 5.0 | | | |

Ref: TEMP-0382 Rev.41 FRONT / RECTO 1300090305

* For research use only in Canada and USA

ABX Minotrol Retic



LOT RX 435

CONTROL

(Exp.) 2022-07-05
(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITS UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|--------------|--|---|-------------------------|--------------|--|---|-------------------------|--------------|--|---|-------------------------|
| | | CONTROL | | | TOLERANCES TOLERANCE | CONTROL | | | TOLERANCES TOLERANCE | CONTROL | | | TOLERANCES TOLERANCE |
| | | PENTRA | | 1 | | PENTRA | | 2 | | PENTRA | | 3 | |
| | | XLR | | | | XLR | | | | XLR | | | |
| Méthode Manuelle / Manual RETIC | % | 1.3 | | | ± 0.7 | 5.2 | | | ± 2.1 | 10.1 | | | ± 3.0 |
| RETIC | % | 1.58 | | | ± 0.79 | 5.20 | | | ± 1.30 | 10.76 | | | ± 2.15 |
| RETIC | 10 ⁶ /µl | 0.053 | | | ± 0.027 | 0.173 | | | ± 0.043 | 0.353 | | | ± 0.071 |
| MFI MRI * | % | 28.1 | | | ± 8.0 | 43.4 | | | ± 10.0 | 46.5 | | | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.37 | | | ± 0.30 | 3.33 | | | ± 0.30 | 3.28 | | | ± 0.30 |
| PIC RANGE | - | 8 | | | ± 4 | 8 | | | ± 4 | 8 | | | ± 4 |

Ref: TEMP-0382 Rev.41 BACK / VERSO 1300090305

* For research use only in Canada and USA