

ABX Minotrol Retic



LOT

RX 090

CONTROL

(Exp.) **2010-11-05**
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		120 RETIC	DX120	120 RETIC	DX120	120 RETIC	DX120						
Méthode Manuelle / Manual													
RETIC	%	1.3	1.3		± 0.7	6.0	6.0		± 2.4	11.2	11.2		± 3.4
RETIC	%	1.22	1.31		± 0.66	5.06	5.10		± 1.28	10.30	10.50		± 2.10
RETIC	10 ⁹ /µl	0.041	0.045		± 0.023	0.175	0.177		± 0.044	0.353	0.364		± 0.073
MFI MRI *		24.9	25.0		± 8.0	38.4	38.6		± 10.0	41.7	41.4		± 10.0
GR RBC	10 ⁹ /µl	3.40	3.40		± 0.30	3.45	3.47		± 0.30	3.43	3.47		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

* For research use only in Canada and USA

9930034-C
Ref: FORM-0382 Rev.06