

ABX Minotrol Retic



HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

RX 012

CONTROL

(Exp.) **2012-03-05**
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		120 RETIC	DX120		120 RETIC	DX120		120 RETIC	DX120				
Méthode Manuelle / Manual RETIC	%	1.4	1.4		± 0.7	6.3	6.3		± 2.5	10.8	10.8		± 3.2
RETIC	%	1.30	1.25		± 0.65	5.33	5.38		± 1.35	10.90	11.00		± 2.20
RETIC	10 ⁶ /µl	0.045	0.044		± 0.023	0.184	0.188		± 0.047	0.374	0.385		± 0.077
MFI MRI *		22.0	21.8		± 8.0	40.4	42.7		± 10.0	43.8	45.4		± 10.0
GR RBC	10 ⁶ /µl	3.43	3.49		± 0.30	3.45	3.49		± 0.30	3.43	3.50		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C

Ref: FORM-0382 Rev.07

* For research use only in Canada and USA