

ABX Minotrol Retic



LOT RX 032
Rev 1

CONTROL

(Exp.) 2012-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		120 RETIC	DX120	120 RETIC	DX120	120 RETIC	DX120						
Méthode Manuelle / Manual													
RETIC	%	1.3	1.3		± 0.7	5.8	5.8		± 2.3	10.9	10.9		± 3.3
RETIC	%	1.68	1.86		± 0.93	5.46	5.58		± 1.40	10.70	10.90		± 2.20
RETIC	10 ⁹ /µl	0.057	0.064		± 0.032	0.187	0.193		± 0.048	0.372	0.380		± 0.076
MFI MRI *		32.1	32.6		± 8.0	44.7	44.7		± 10.0	48.8	49.1		± 10.0
GR RBC	10 ⁶ /µl	3.40	3.43		± 0.30	3.42	3.45		± 0.30	3.48	3.49		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C

Ref: FORM-0382 Rev.07

* For research use only in Canada and USA