

# ABX Minotrol Retic



**LOT**

RX 111

**CONTROL**

(Exp.) 2012-01-05  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA	PENTRA		PENTRA	PENTRA	PENTRA					
		120 RETIC	DX120	120 RETIC	DX120	120 RETIC	DX120						
Méthode Manuelle / Manual RETIC	%	1.3	1.3		± 0.7	6.4	6.4		± 2.6	11.4	11.4		± 3.4
RETIC	%	1.26	1.54		± 0.77	5.22	5.47		± 1.37	10.40	10.60		± 2.20
RETIC	10 <sup>9</sup> /µl	0.044	0.054		± 0.027	0.180	0.190		± 0.048	0.358	0.368		± 0.074
MFI MRI *		26.4	24.9		± 8.0	41.5	41.9		± 10.0	43.2	45.9		± 10.0
GR RBC	10 <sup>9</sup> /µl	3.47	3.50		± 0.30	3.44	3.47		± 0.30	3.44	3.47		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C  
Ref: FORM-0382 Rev.07

\* For research use only in Canada and USA