

# ABX Minotrol Retic



**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

RX 071

**CONTROL**

(Exp.) **2011-09-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA	PENTRA		PENTRA	PENTRA	PENTRA					
		120 RETIC	DX120	120 RETIC	DX120	120 RETIC	DX120						
Méthode Manuelle / Manual RETIC	%	1.4	1.4		± 0.7	6.0	6.0		± 2.5	11.0	11.0		± 3.5
RETIC	%	1.38	1.43		± 0.72	5.28	5.36		± 1.34	11.00	11.00		± 2.20
RETIC	10 <sup>9</sup> /µl	0.047	0.048		± 0.024	0.181	0.181		± 0.045	0.385	0.380		± 0.077
MFI MRI *		29.7	32.4		± 8.0	42.4	43.9		± 10.0	45.8	47.2		± 10.0
GR RBC	10 <sup>9</sup> /µl	3.43	3.37		± 0.30	3.43	3.37		± 0.30	3.50	3.45		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C

Ref: FORM-0382 Rev.07

\* For research use only in Canada and USA