

# ABX Minotrol Retic



**LOT** RX 056  
Rev 1

**CONTROL**

(Exp.) 2016-07-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX       |        |          |                         |           |        |          |                         |           |       |          |                         |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
|                                       |                     | CONTROL 1 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 2 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 3 |       |          | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    |       |          |                         |
|                                       |                     | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120 | DX Nexus |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | 1.6       | 1.6    | 1.6      | ± 0.8                   | 5.5       | 5.5    | 5.5      | ± 2.2                   | 10.2      | 10.2  | 10.2     | ± 3.1                   |
| RETIC                                 | %                   | 2.10      | 2.11   | 2.11     | ± 1.06                  | 5.66      | 5.66   | 5.66     | ± 1.42                  | 11.70     | 11.60 | 11.60    | ± 2.40                  |
| RETIC                                 | 10 <sup>6</sup> /µl | 0.073     | 0.072  | 0.072    | ± 0.037                 | 0.195     | 0.195  | 0.195    | ± 0.049                 | 0.401     | 0.394 | 0.394    | ± 0.082                 |
| MFI MRI *                             |                     | 30.6      | 31.0   | 31.0     | ± 8.0                   | 44.7      | 44.8   | 44.8     | ± 10.0                  | 49.2      | 50.6  | 50.6     | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | 3.47      | 3.43   | 3.43     | ± 0.30                  | 3.45      | 3.45   | 3.45     | ± 0.30                  | 3.43      | 3.40  | 3.40     | ± 0.30                  |
| PIC RANGE                             |                     | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0  | 27.0     | ± 5.0                   |

Ref: TEMP-0382 Rev.34 FRONT / RECTO 9930083 -A

\* For research use only in Canada and USA

# ABX Minotrol Retic



**LOT** RX 056  
Rev 1

**CONTROL**

(Exp.) 2016-07-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX          |        |        |                         |              |        |        |                         |              |        |  |                         |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
|                                       |                     | CONTROL 1    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 2    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 3    |        |  | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA |  |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | <b>1.6</b>   |        |        | ± 0.8                   | <b>5.5</b>   |        |        | ± 2.2                   | <b>10.2</b>  |        |  | ± 3.1                   |
| RETIC                                 | %                   | <b>2.00</b>  |        |        | ± 1.00                  | <b>5.47</b>  |        |        | ± 1.37                  | <b>11.30</b> |        |  | ± 2.26                  |
| RETIC                                 | 10 <sup>6</sup> /µl | <b>0.069</b> |        |        | ± 0.035                 | <b>0.187</b> |        |        | ± 0.047                 | <b>0.385</b> |        |  | ± 0.077                 |
| MFI MRI *                             |                     | <b>28.9</b>  |        |        | ± 8.0                   | <b>43.1</b>  |        |        | ± 10.0                  | <b>49.1</b>  |        |  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | <b>3.43</b>  |        |        | ± 0.30                  | <b>3.41</b>  |        |        | ± 0.30                  | <b>3.41</b>  |        |  | ± 0.30                  |
| PIC RANGE                             |                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |  | ± 4                     |

Ref: TEMP-0382 Rev.34 BACK / VERSO 9930083-A

\* For research use only in Canada and USA