

# ABX Minotrol Retic



**LOT** RX 095  
Rev 1

**CONTROL**

(Exp.) 2015-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX       |           |        |                         |           |        |          |                         |           |        |       |                         |
|---------------------------------------|---------------------|-----------|-----------|--------|-------------------------|-----------|--------|----------|-------------------------|-----------|--------|-------|-------------------------|
|                                       |                     | CONTROL 1 |           |        | TOLERANCES<br>TOLERANCE | CONTROL 2 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 3 |        |       | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA    | PENTRA    | PENTRA |                         | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    | PENTRA |       |                         |
| 120 RETIC                             | DX120               | DX Nexus  | 120 RETIC | DX120  | DX Nexus                | 120 RETIC | DX120  | DX Nexus |                         |           |        |       |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | 1.5       | 1.5       | 1.5    | ± 0.8                   | 5.7       | 5.7    | 5.7      | ± 2.3                   | 10.3      | 10.3   | 10.3  | ± 3.1                   |
| RETIC                                 | %                   | 2.43      | 2.44      | 2.44   | ± 1.22                  | 5.65      | 5.87   | 5.87     | ± 1.47                  | 11.70     | 11.80  | 11.80 | ± 2.40                  |
| RETIC                                 | 10 <sup>6</sup> /µl | 0.083     | 0.084     | 0.084  | ± 0.042                 | 0.194     | 0.203  | 0.203    | ± 0.051                 | 0.401     | 0.405  | 0.405 | ± 0.081                 |
| MFI MRI *                             |                     | 33.5      | 33.5      | 33.5   | ± 8.0                   | 45.1      | 44.5   | 44.5     | ± 10.0                  | 49.9      | 48.9   | 48.9  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | 3.42      | 3.45      | 3.45   | ± 0.30                  | 3.44      | 3.45   | 3.45     | ± 0.30                  | 3.43      | 3.43   | 3.43  | ± 0.30                  |
| PIC RANGE                             |                     | 27.0      | 27.0      | 27.0   | ± 5.0                   | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0   | 27.0  | ± 5.0                   |

Ref: TEMP-0382 Rev.34 FRONT / RECTO 9930083 -A

\* For research use only in Canada and USA

# ABX Minotrol Retic



**LOT** RX 095  
Rev 1

**CONTROL**

(Exp.) 2015-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX          |        |        |                         |              |        |        |                         |              |        |  |                         |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
|                                       |                     | CONTROL 1    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 2    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 3    |        |  | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA |  |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | <b>1.5</b>   |        |        | ± 0.8                   | <b>5.7</b>   |        |        | ± 2.3                   | <b>10.3</b>  |        |  | ± 3.1                   |
| RETIC                                 | %                   | <b>2.44</b>  |        |        | ± 1.22                  | <b>5.67</b>  |        |        | ± 1.42                  | <b>11.50</b> |        |  | ± 2.30                  |
| RETIC                                 | 10 <sup>6</sup> /µl | <b>0.083</b> |        |        | ± 0.042                 | <b>0.193</b> |        |        | ± 0.048                 | <b>0.390</b> |        |  | ± 0.078                 |
| MFI MRI *                             |                     | <b>35.9</b>  |        |        | ± 8.0                   | <b>46.4</b>  |        |        | ± 10.0                  | <b>49.8</b>  |        |  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | <b>3.41</b>  |        |        | ± 0.30                  | <b>3.40</b>  |        |        | ± 0.30                  | <b>3.39</b>  |        |  | ± 0.30                  |
| PIC RANGE                             |                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |  | ± 4                     |

Ref: TEMP-0382 Rev.34 BACK / VERSO 9930083-A

\* For research use only in Canada and USA