

# ABX Minotrol Retic



**LOT** RX 419  
Rev 1

**CONTROL**

(Exp.) 2019-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX       |        |          |                         |           |        |          |                         |           |       |          |                         |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
|                                       |                     | CONTROL 1 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 2 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 3 |       |          | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    |       |          |                         |
|                                       |                     | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120 | DX Nexus |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | 1.5       | 1.5    | 1.5      | ± 0.8                   | 4.8       | 4.8    | 4.8      | ± 2.1                   | 10.5      | 10.5  | 10.5     | ± 3.2                   |
| RETIC                                 | %                   | 1.44      | 1.64   | 1.64     | ± 0.82                  | 4.72      | 5.20   | 5.20     | ± 1.30                  | 10.00     | 10.90 | 10.90    | ± 2.20                  |
| RETIC                                 | 10 <sup>6</sup> /µl | 0.049     | 0.056  | 0.056    | ± 0.028                 | 0.160     | 0.176  | 0.176    | ± 0.044                 | 0.343     | 0.375 | 0.375    | ± 0.075                 |
| MFI MRI *                             |                     | 30.0      | 27.9   | 27.9     | ± 8.0                   | 43.0      | 41.3   | 41.3     | ± 10.0                  | 46.2      | 44.0  | 44.0     | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | 3.42      | 3.42   | 3.42     | ± 0.30                  | 3.38      | 3.39   | 3.39     | ± 0.30                  | 3.43      | 3.44  | 3.44     | ± 0.30                  |
| PIC RANGE                             |                     | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0  | 27.0     | ± 5.0                   |

\* For research use only in Canada and USA

Ref: TEMP-0382 Rev.36 FRONT / RECTO 1219930083-A

# ABX Minotrol Retic



**LOT** RX 419  
Rev 1

**CONTROL**

(Exp.) 2019-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX          |        |        |                         |              |        |        |                         |              |        |  |                         |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
|                                       |                     | CONTROL 1    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 2    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 3    |        |  | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA |  |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | <b>1.5</b>   |        |        | ± 0.8                   | <b>4.8</b>   |        |        | ± 2.1                   | <b>10.5</b>  |        |  | ± 3.2                   |
| RETIC                                 | %                   | <b>1.62</b>  |        |        | ± 0.81                  | <b>5.14</b>  |        |        | ± 1.29                  | <b>10.65</b> |        |  | ± 2.15                  |
| RETIC                                 | 10 <sup>6</sup> /µl | <b>0.055</b> |        |        | ± 0.028                 | <b>0.172</b> |        |        | ± 0.043                 | <b>0.362</b> |        |  | ± 0.072                 |
| MFI MRI *                             |                     | <b>24.9</b>  |        |        | ± 8.0                   | <b>39.6</b>  |        |        | ± 10.0                  | <b>43.8</b>  |        |  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | <b>3.40</b>  |        |        | ± 0.30                  | <b>3.35</b>  |        |        | ± 0.30                  | <b>3.40</b>  |        |  | ± 0.30                  |
| PIC RANGE                             |                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |  | ± 4                     |

Ref: TEMP-0382 Rev.36 BACK / VERSO 1219930083-A

\* For research use only in Canada and USA