

# ABX Minotrol Retic



**LOT** RX 413  
Rev 1

**CONTROL**

(Exp.) 2018-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX       |        |          |                         |           |        |          |                         |           |       |          |                         |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
|                                       |                     | CONTROL 1 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 2 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 3 |       |          | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    |       |          |                         |
|                                       |                     | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120 | DX Nexus |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | 1.3       | 1.3    | 1.3      | ± 0.7                   | 4.8       | 4.8    | 4.8      | ± 2.0                   | 9.8       | 9.8   | 9.8      | ± 3.0                   |
| RETIC                                 | %                   | 1.82      | 1.86   | 1.86     | ± 0.93                  | 5.11      | 5.23   | 5.23     | ± 1.31                  | 10.90     | 11.40 | 11.40    | ± 2.30                  |
| RETIC                                 | 10 <sup>6</sup> /µl | 0.064     | 0.063  | 0.063    | ± 0.032                 | 0.182     | 0.179  | 0.179    | ± 0.046                 | 0.379     | 0.380 | 0.380    | ± 0.076                 |
| MFI MRI *                             |                     | 29.9      | 27.3   | 27.3     | ± 8.0                   | 48.1      | 45.4   | 45.4     | ± 10.0                  | 51.0      | 49.5  | 49.5     | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | 3.54      | 3.41   | 3.41     | ± 0.30                  | 3.57      | 3.43   | 3.43     | ± 0.30                  | 3.48      | 3.33  | 3.33     | ± 0.30                  |
| PIC RANGE                             |                     | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0  | 27.0     | ± 5.0                   |

Ref: TEMP-0382 Rev.36 FRONT / RECTO 1219930083-A

\* For research use only in Canada and USA

# ABX Minotrol Retic



**LOT** RX 413  
Rev 1

**CONTROL**

(Exp.) 2018-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX          |        |        |                         |              |        |        |                         |              |        |  |                         |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
|                                       |                     | CONTROL 1    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 2    |        |        | TOLERANCES<br>TOLERANCE | CONTROL 3    |        |  | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA |  |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | <b>1.3</b>   |        |        | ± 0.7                   | <b>4.8</b>   |        |        | ± 2.0                   | <b>9.8</b>   |        |  | ± 3.0                   |
| RETIC                                 | %                   | <b>1.78</b>  |        |        | ± 0.89                  | <b>5.04</b>  |        |        | ± 1.26                  | <b>11.00</b> |        |  | ± 2.20                  |
| RETIC                                 | 10 <sup>6</sup> /µl | <b>0.061</b> |        |        | ± 0.031                 | <b>0.174</b> |        |        | ± 0.044                 | <b>0.366</b> |        |  | ± 0.073                 |
| MFI MRI *                             |                     | <b>27.7</b>  |        |        | ± 8.0                   | <b>45.9</b>  |        |        | ± 10.0                  | <b>51.3</b>  |        |  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | <b>3.42</b>  |        |        | ± 0.30                  | <b>3.45</b>  |        |        | ± 0.30                  | <b>3.33</b>  |        |  | ± 0.30                  |
| PIC RANGE                             |                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |  | ± 4                     |

Ref: TEMP-0382 Rev.36 BACK / VERSO 1219930083-A

\* For research use only in Canada and USA