

ABX Minotrol Retic



LOT RX 404
Rev 2

CONTROL

(Exp.) 2017-05-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | | | |
| | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | |
| Méthode Manuelle / Manual RETIC | % | 1.7 | 1.7 | 1.7 | ± 0.9 | 5.5 | 5.5 | 5.5 | ± 2.2 | 12.0 | 12.0 | 12.0 | ± 4.2 |
| RETIC | % | 1.95 | 2.13 | 2.13 | ± 1.07 | 5.38 | 5.48 | 5.48 | ± 1.37 | 11.00 | 11.00 | 11.00 | ± 2.20 |
| RETIC | 10 ⁶ /µl | 0.067 | 0.073 | 0.073 | ± 0.037 | 0.182 | 0.184 | 0.184 | ± 0.046 | 0.372 | 0.367 | 0.367 | ± 0.074 |
| MFI MRI * | | 27.2 | 27.6 | 27.6 | ± 8.0 | 43.4 | 44.9 | 44.9 | ± 10.0 | 46.3 | 49.4 | 49.4 | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.45 | 3.44 | 3.44 | ± 0.30 | 3.39 | 3.36 | 3.36 | ± 0.30 | 3.38 | 3.34 | 3.34 | ± 0.30 |
| PIC RANGE | | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 |

Ref: TEMP-0382 Rev.36 FRONT / RECTO 1219930083-A

* For research use only in Canada and USA

ABX Minotrol Retic



LOT RX 404
Rev 2

CONTROL

(Exp.) 2017-05-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | | |
| Méthode Manuelle / Manual RETIC | % | 1.7 | | | ± 0.9 | 5.5 | | | ± 2.2 | 12.0 | | | ± 4.2 |
| RETIC | % | 1.89 | | | ± 0.95 | 5.28 | | | ± 1.32 | 11.00 | | | ± 2.20 |
| RETIC | 10 ⁶ /µl | 0.067 | | | ± 0.034 | 0.181 | | | ± 0.045 | 0.375 | | | ± 0.075 |
| MFI MRI * | | 28.8 | | | ± 8.0 | 45.3 | | | ± 10.0 | 48.7 | | | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.52 | | | ± 0.30 | 3.42 | | | ± 0.30 | 3.41 | | | ± 0.30 |
| PIC RANGE | | 8 | | | ± 4 | 8 | | | ± 4 | 8 | | | ± 4 |

Ref: TEMP-0382 Rev.36 BACK / VERSO 1219930083-A

* For research use only in Canada and USA