

# ABX Minotrol Retic



**LOT** RX 401  
Rev 1

**CONTROL**

(Exp.) 2016-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX       |        |          |                         |           |        |          |                         |           |       |          |                         |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
|                                       |                     | CONTROL 1 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 2 |        |          | TOLERANCES<br>TOLERANCE | CONTROL 3 |       |          | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    | PENTRA | PENTRA   |                         | PENTRA    |       |          |                         |
|                                       |                     | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120  | DX Nexus |                         | 120 RETIC | DX120 | DX Nexus |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | 1.8       | 1.8    | 1.8      | ± 1.0                   | 6.1       | 6.1    | 6.1      | ± 2.4                   | 11.0      | 11.0  | 11.0     | ± 3.3                   |
| RETIC                                 | %                   | 2.14      | 2.15   | 2.15     | ± 1.08                  | 5.69      | 5.78   | 5.78     | ± 1.45                  | 11.50     | 11.50 | 11.50    | ± 2.30                  |
| RETIC                                 | 10 <sup>6</sup> /µl | 0.075     | 0.075  | 0.075    | ± 0.038                 | 0.197     | 0.200  | 0.200    | ± 0.050                 | 0.397     | 0.394 | 0.394    | ± 0.079                 |
| MFI MRI *                             |                     | 27.6      | 27.3   | 27.3     | ± 8.0                   | 47.7      | 46.7   | 46.7     | ± 10.0                  | 49.9      | 48.9  | 48.9     | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | 3.51      | 3.50   | 3.50     | ± 0.30                  | 3.47      | 3.46   | 3.46     | ± 0.30                  | 3.45      | 3.43  | 3.43     | ± 0.30                  |
| PIC RANGE                             |                     | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0   | 27.0     | ± 5.0                   | 27.0      | 27.0  | 27.0     | ± 5.0                   |

Ref: TEMP-0382 Rev.34 FRONT / RECTO 9930083 -A

\* For research use only in Canada and USA

# ABX Minotrol Retic



**LOT** RX 401  
Rev 1

**CONTROL**

(Exp.) 2016-11-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS              | UNITES<br>UNITS     | ABX          |        |        |                         |              |        |        |                         |              |        |  |                         |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--------|--|-------------------------|
|                                       |                     | CONTROL      |        |        | TOLERANCES<br>TOLERANCE | CONTROL      |        |        | TOLERANCES<br>TOLERANCE | CONTROL      |        |  | TOLERANCES<br>TOLERANCE |
|                                       |                     | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA | PENTRA |                         | PENTRA       | PENTRA |  |                         |
| Méthode<br>Manuelle / Manual<br>RETIC | %                   | <b>1.8</b>   |        |        | ± 1.0                   | <b>6.1</b>   |        |        | ± 2.4                   | <b>11.0</b>  |        |  | ± 3.3                   |
| RETIC                                 | %                   | <b>2.14</b>  |        |        | ± 1.07                  | <b>5.61</b>  |        |        | ± 1.40                  | <b>11.30</b> |        |  | ± 2.26                  |
| RETIC                                 | 10 <sup>6</sup> /µl | <b>0.075</b> |        |        | ± 0.038                 | <b>0.195</b> |        |        | ± 0.049                 | <b>0.390</b> |        |  | ± 0.078                 |
| MFI MRI *                             |                     | <b>26.7</b>  |        |        | ± 8.0                   | <b>45.6</b>  |        |        | ± 10.0                  | <b>47.6</b>  |        |  | ± 12.0                  |
| GR RBC                                | 10 <sup>6</sup> /µl | <b>3.52</b>  |        |        | ± 0.30                  | <b>3.47</b>  |        |        | ± 0.30                  | <b>3.45</b>  |        |  | ± 0.30                  |
| PIC RANGE                             |                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |        | ± 4                     | <b>8</b>     |        |  | ± 4                     |

Ref: TEMP-0382 Rev.34 BACK / VERSO 9930083-A

\* For research use only in Canada and USA