

ABX Minotrol Retic



LOT RX 407
Rev 1

CONTROL

(Exp.) 2017-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.4			± 0.7	5.9			± 2.4	11.2			± 3.4
RETIC	%	2.11			± 1.06	5.58			± 1.40	11.20			± 2.30
RETIC	10 ⁶ /mm ³	0.073			± 0.037	0.191			± 0.048	0.391			± 0.078
	10 ⁹ /l	73			± 37	191			± 48	391			± 78
MFI MRI *		27.8			± 8.0	39.3			± 10.0	42.3			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.45			± 0.30	3.43			± 0.30	3.49			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA