

ABX Minotrol Retic



LOT RX 418

CONTROL

(Exp.) 2019-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.5			± 0.8	5.6			± 2.2	11.0			± 3.3
RETIC	%	1.68			± 0.84	5.19			± 1.30	11.10			± 2.30
RETIC	10 ⁶ /mm ³	0.059			± 0.030	0.182			± 0.046	0.391			± 0.078
	10 ⁹ /l	59			± 30	182			± 46	391			± 78
MFI MRI *		29.7			± 8.0	39.4			± 10.0	43.2			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.52			± 0.30	3.51			± 0.30	3.52			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref. TEMP-0382 Rev.36

* For research use only in Canada and USA