

ABX Minotrol Retic



LOT RX 417
Rev 1

CONTROL

(Exp.) 2019-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	5.4			± 2.2	10.9			± 3.3
RETIC	%	1.84			± 0.92	5.76			± 1.44	11.15			± 2.23
RETIC	10 ⁶ /mm ³	0.065			± 0.033	0.201			± 0.051	0.389			± 0.078
	10 ⁹ /l	65			± 33	201			± 51	389			± 78
MFI MRI *		32.8			± 8.0	47.2			± 10.0	49.0			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.52			± 0.30	3.49			± 0.30	3.49			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

* For research use only in Canada and USA

1300032515-A
Ref: TEMP-0382 Rev.36