

ABX Minotrol Retic



LOT RX 426
Rev 1

CONTROL

(Exp.) 2021-01-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.4			± 0.8	5.8			± 2.3	11.6			± 3.8
RETIC	%	1.77			± 0.89	5.48			± 1.37	10.88			± 2.18
RETIC	10 ⁶ /mm ³	0.062			± 0.031	0.193			± 0.048	0.384			± 0.077
	10 ⁹ /l	62			± 31	193			± 48	384			± 77
MFI MRI *	%	30.2			± 8.0	39.0			± 10.0	41.0			± 12.0
	-	0.302			± 0.080	0.390			± 0.100	0.410			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.50			± 0.30	3.52			± 0.30	3.53			± 0.30
PIC RANGE	-	13			± 4	13			± 4	13			± 4

1300080026
Ref: TEMP-0382 Rev.38

* For research use only in Canada and USA