

ABX Minotrol Retic



LOT RX 416
Rev 1

CONTROL

(Exp.) 2019-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.4			± 0.7	4.7			± 1.9	11.8			± 3.5
RETIC	%	1.61			± 0.81	5.47			± 1.37	11.30			± 2.30
RETIC	10 ⁶ /mm ³	0.057			± 0.029	0.195			± 0.049	0.397			± 0.079
	10 ⁹ /l	57			± 29	195			± 49	397			± 79
MFI MRI *		25.9			± 8.0	40.4			± 10.0	43.3			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.53			± 0.30	3.57			± 0.30	3.51			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA