

ABX Minotrol Retic



LOT RX 406
Rev 1

CONTROL

(Exp.) 2017-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.5			± 1.0	5.7			± 2.3	11.6			± 3.5
RETIC	%	1.90			± 0.95	5.63			± 1.41	12.00			± 2.40
RETIC	10 ⁶ /mm ³	0.068			± 0.034	0.199			± 0.050	0.424			± 0.085
	10 ⁹ /l	68			± 34	199			± 50	424			± 85
MFI MRI *		25.8			± 8.0	43.3			± 10.0	49.0			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.56			± 0.30	3.54			± 0.30	3.53			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

* For research use only in Canada and USA

1300032515-A
Ref: TEMP-0382 Rev.36