

ABX Minotrol Retic



LOT RX 415
Rev 3

CONTROL

(Exp.) 2019-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
H2500			H2500			H2500							
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	4.9			± 2.0	10.0			± 3.0
RETIC	%	1.41			± 0.71	5.02			± 1.26	11.00			± 2.20
RETIC	10 ⁶ /mm ³	0.050			± 0.025	0.180			± 0.045	0.394			± 0.079
	10 ⁹ /l	50			± 25	180			± 45	394			± 79
MFI MRI *		26.9			± 8.0	41.2			± 10.0	44.1			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.55			± 0.30	3.58			± 0.30	3.58			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA