

ABX Minotrol Retic



LOT RX 405
Rev 1

CONTROL

(Exp.) 2017-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.5			± 0.8	6.0			± 2.4	11.7			± 3.5
RETIC	%	2.05			± 1.03	5.80			± 1.45	12.00			± 2.40
RETIC	10 ⁶ /mm ³	0.071			± 0.036	0.202			± 0.051	0.413			± 0.083
	10 ⁹ /l	71			± 36	202			± 51	413			± 83
MFI MRI *		25.0			± 8.0	44.6			± 10.0	49.4			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.47			± 0.30	3.48			± 0.30	3.44			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

* For research use only in Canada and USA

1300032515-A
Ref: TEMP-0382 Rev.36