

ABX Minotrol Retic



LOT RX 424
Rev 1

CONTROL

(Exp.) 2020-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.3			± 1.0	4.9			± 2.4	10.7			± 3.2
RETIC	%	1.58			± 0.79	4.92			± 1.23	10.89			± 2.18
RETIC	10 ⁶ /mm ³	0.056			± 0.028	0.173			± 0.043	0.374			± 0.075
	10 ⁹ /l	56			± 28	173			± 43	374			± 75
MFI MRI *	%	26.9			± 8.0	41.3			± 10.0	44.0			± 12.0
	-	0.269			± 0.080	0.413			± 0.100	0.440			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.53			± 0.30	3.52			± 0.30	3.43			± 0.30
PIC RANGE	-	13			± 4	13			± 4	13			± 4

1300080026
Ref: TEMP-0382 Rev.37

* For research use only in Canada and USA