

# ABX Minotrol Retic



**LOT** RX 421  
Rev 1

**CONTROL**

(Exp.) 2020-03-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.7			± 0.9	7.1			± 2.8	11.3			± 3.4
RETIC	%	1.94			± 0.97	5.64			± 1.41	11.45			± 2.29
RETIC	10 <sup>6</sup> /mm <sup>3</sup>	0.068			± 0.034	0.196			± 0.049	0.403			± 0.081
	10 <sup>9</sup> /l	68			± 34	196			± 49	403			± 81
MFI MRI *		32.9			± 8.0	40.8			± 10.0	44.1			± 12.0
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	3.48			± 0.30	3.48			± 0.30	3.52			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A  
Ref: TEMP-0382 Rev.36

\* For research use only in Canada and USA