

ABX Minotrol Retic



LOT RX 411
Rev 1

CONTROL

(Exp.) 2018-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.4			± 0.8	6.0			± 2.4	10.8			± 3.2
RETIC	%	2.10			± 1.05	5.54			± 1.39	11.60			± 2.32
RETIC	10 ⁶ /mm ³	0.073			± 0.037	0.193			± 0.048	0.404			± 0.081
	10 ⁹ /l	73			± 37	193			± 48	404			± 81
MFI MRI *		34.7			± 8.0	43.9			± 10.0	48.0			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.49			± 0.30	3.48			± 0.30	3.48			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA