

ABX Minotrol Retic



LOT RX 420
Rev 1

CONTROL

(Exp.) 2020-01-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.4			± 0.7	5.3			± 2.1	11.3			± 3.4
RETIC	%	1.87			± 0.94	5.62			± 1.41	11.16			± 2.23
RETIC	10 ⁶ /mm ³	0.066			± 0.033	0.199			± 0.050	0.391			± 0.078
	10 ⁹ /l	66			± 33	199			± 50	391			± 78
MFI MRI *		34.8			± 8.0	42.3			± 10.0	44.8			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.52			± 0.30	3.54			± 0.30	3.50			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

Ref: TEMP-0382 Rev.36 1300032515-A

* For research use only in Canada and USA