

ABX Minotrol Retic



LOT RX 410
Rev 1

CONTROL

(Exp.) 2018-05-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---|--------------|--|---|-------------------------|--------------|--|---|-------------------------|--------------|--|---|-------------------------|
| | | CONTROL | | | TOLERANCES TOLERANCE | CONTROL | | | TOLERANCES TOLERANCE | CONTROL | | | TOLERANCES TOLERANCE |
| | | YUMIZEN | | 1 | | YUMIZEN | | 2 | | YUMIZEN | | 3 | |
| | | H2500 | | | | H2500 | | | | H2500 | | | |
| Méthode Manuelle / Manual RETIC | % | 1.5 | | | ± 0.8 | 6.1 | | | ± 2.4 | 10.5 | | | ± 3.2 |
| RETIC | % | 1.80 | | | ± 0.90 | 5.92 | | | ± 1.48 | 11.30 | | | ± 2.26 |
| RETIC | 10 ⁶ /mm ³ | 0.062 | | | ± 0.031 | 0.205 | | | ± 0.051 | 0.374 | | | ± 0.075 |
| | 10 ⁹ /l | 62 | | | ± 31 | 205 | | | ± 51 | 374 | | | ± 75 |
| MFI MRI * | | 31.7 | | | ± 8.0 | 45.4 | | | ± 10.0 | 48.9 | | | ± 12.0 |
| GR RBC | 10 ⁶ /mm ³ ; 10 ¹² /l | 3.43 | | | ± 0.30 | 3.46 | | | ± 0.30 | 3.31 | | | ± 0.30 |
| PIC RANGE | | 13 | | | ± 4 | 13 | | | ± 4 | 13 | | | ± 4 |

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA