

ABX Minotrol Retic



LOT RX 409
Rev 1

CONTROL

(Exp.) 2018-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.5			± 0.9	5.5			± 2.2	11.4			± 3.4
RETIC	%	1.88			± 0.94	5.37			± 1.34	11.50			± 2.30
RETIC	10 ⁶ /mm ³	0.065			± 0.033	0.192			± 0.048	0.413			± 0.083
	10 ⁹ /l	65			± 33	192			± 48	413			± 83
MFI MRI *		27.1			± 8.0	43.6			± 10.0	50.3			± 12.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.48			± 0.30	3.57			± 0.30	3.59			± 0.30
PIC RANGE		13			± 4	13			± 4	13			± 4

1300032515-A
Ref: TEMP-0382 Rev.36

* For research use only in Canada and USA