

ABX Minotrol Retic



LOT RX 433
Rev 1

CONTROL

(Exp.) 2022-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.7			± 0.9	5.5			± 2.2	11.0			± 3.3
RETIC	%	1.56			± 0.78	5.13			± 1.28	10.40			± 2.08
RETIC	10 ⁶ /mm ³	0.054			± 0.027	0.177			± 0.044	0.362			± 0.072
	10 ⁹ /l	54			± 27	177			± 44	362			± 72
MFI MRI *	%	31.4			± 8.0	43.7			± 10.0	47.1			± 12.0
	-	0.314			± 0.080	0.437			± 0.100	0.471			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.45			± 0.30	3.46			± 0.30	3.48			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

1300080026
Ref: TEMP-0382 Rev.40

* For research use only in Canada and USA