

ABX Minotrol Retic



LOT RX 435
Rev 1

CONTROL

(Exp.) 2022-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	5.2			± 2.1	10.1			± 3.0
RETIC	%	1.60			± 0.80	5.27			± 1.32	11.08			± 2.22
RETIC	10 ⁶ /mm ³	0.057			± 0.029	0.184			± 0.046	0.383			± 0.077
	10 ⁹ /l	57			± 29	184			± 46	383			± 77
MFI MRI *	%	29.7			± 8.0	46.5			± 10.0	49.0			± 12.0
	-	0.297			± 0.080	0.465			± 0.100	0.490			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.54			± 0.30	3.50			± 0.30	3.46			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

1300102875
Ref: TEMP-0382 Rev.41

* For research use only in Canada and USA