


LOT PX 403
Rev 2


CONTROL

 (Exp.) 2017-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio																		
			CONTROL					L	CONTROL					N	CONTROL					H	TOLERANCES TOLERANCE
			PENTRA						PENTRA						PENTRA						
			60 60C+ ES60	80 XL80	MS60	XLR	MSCR		60 60C+ ES60	80 XL80	MS60	XLR	MSCR		60 60C+ ES60	80 XL80	MS60	XLR	MSCR		
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.5	2.4	2.5	2.5	± 0.4	7.5	7.5	7.4	7.5	7.6	± 1.0	18.6	18.2	18.2	18.2	18.7	± 2.2	
GR	RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.38	2.34	2.37	2.34	2.34	± 0.16	4.60	4.62	4.59	4.62	4.57	± 0.20	5.10	5.13	5.10	5.13	5.10	± 0.25	
HB	HGB	g/dl	6.9	6.8	6.8	6.8	6.8	± 0.4	13.8	13.8	13.7	13.8	13.7	± 0.5	16.2	16.3	16.3	16.3	16.2	± 0.6	
		g/l	69	68	68	68	68	± 4	138	138	137	138	137	± 5	162	163	163	163	162	± 6	
		mmol/l	4.28	4.22	4.22	4.22	4.22	± 0.25	8.57	8.57	8.51	8.57	8.51	± 0.31	10.06	10.12	10.12	10.12	10.06	± 0.37	
HT	HCT	%	19.0	19.4	19.0	19.4	18.8	± 1.5	37.3	37.4	37.2	37.4	37.0	± 2.0	44.4	44.6	43.9	44.6	43.9	± 2.5	
		l/l	0.190	0.194	0.190	0.194	0.188	± 0.015	0.373	0.374	0.372	0.374	0.370	± 0.020	0.444	0.446	0.439	0.446	0.439	± 0.025	
VGM	MCV	µm ³ ; fl	80	83	80	83	80.5	± 5	81	81	81	81	81.0	± 5	87	87	86	87	86.0	± 5	
TGMH	MCH	pg	29.0	29.1	28.7	29.1	29.1	± 2.0	30.0	29.9	29.8	29.9	30.0	± 2.0	31.8	31.8	32.0	31.8	31.8	± 2.5	
		fmol	1.80	1.80	1.78	1.80	1.80	± 0.12	1.86	1.85	1.85	1.85	1.86	± 0.12	1.97	1.97	1.98	1.97	1.97	± 0.16	
CCMH	MCHC	g/dl	36.2	35.0	35.9	35.0	36.1	± 3.0	37.0	36.9	36.8	36.9	37.0	± 3.0	36.5	36.5	37.2	36.5	36.9	± 3.0	
		g/l	362	350	359	350	361	± 30	370	369	368	369	370	± 30	365	365	372	365	369	± 30	
		mmol/l	22.50	21.74	22.27	21.74	22.42	± 1.86	23.00	22.90	22.88	22.90	22.98	± 1.86	22.67	22.68	23.08	22.68	22.94	± 1.86	
IDR	RDW	%	12.5	12.8	12.5	12.8	11.8	± 4.0	13.1	13.2	12.5	13.2	12.5	± 4.0	12.7	13.2	12.0	13.2	12.3	± 4.0	
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	72	67	70	67	68	± 20	242	237	245	237	232	± 30	480	480	495	480	465	± 50	
VPM	MPV	µm ³ ; fl	9.6	10.0	9.7	10.0	9.3	± 2.0	9.5	9.8	9.6	9.8	9.0	± 2.0	9.5	9.8	9.6	9.8	9.1	± 2.0	
NEUT	#		1.42	1.47	1.39	1.47	1.50	± 0.35	4.13	4.16	4.09	4.16	4.20	± 0.90	12.91	12.92	12.78	12.92	13.46	± 1.90	
		%	59.0	58.9	57.8	58.9	59.8	± 10.0	55.1	55.5	55.3	55.5	55.3	± 10.0	69.4	71.0	70.2	71.0	72.0	± 10.0	
LYMPHO	#		0.69	0.73	0.71	0.73	0.73	± 0.33	2.54	2.50	2.53	2.50	2.60	± 0.70	2.88	2.69	2.80	2.69	2.90	± 1.50	
		%	28.8	29.0	29.5	29.0	29.3	± 12.0	33.8	33.3	34.2	33.3	34.2	± 8.0	15.5	14.8	15.4	14.8	15.5	± 8.0	
MONO	#		0.07	0.08	0.09	0.08	0.08	± 0.07	0.32	0.31	0.26	0.31	0.30	± 0.26	0.95	0.87	0.87	0.87	0.75	± 0.75	
		%	3.0	3.1	3.8	3.1	3.0	± 3.0	4.2	4.1	3.5	4.1	4.0	± 3.5	5.1	4.8	4.8	4.8	4.0	± 4.0	
EOS	#		0.14	0.14	0.13	0.14	0.11	± 0.11	0.26	0.27	0.26	0.27	0.23	± 0.23	1.04	0.93	0.95	0.93	0.80	± 0.80	
		%	5.9	5.7	5.5	5.7	4.4	± 4.4	3.4	3.6	3.5	3.6	3.0	± 3.0	5.6	5.1	5.2	5.1	4.3	± 4.3	
BASO	#		0.08	0.08	0.08	0.08	0.09	± 0.08	0.26	0.26	0.26	0.26	0.27	± 0.26	0.82	0.78	0.80	0.78	0.79	± 0.78	
		%	3.3	3.3	3.4	3.3	3.5	± 3.3	3.5	3.5	3.5	3.5	3.5	± 3.5	4.4	4.3	4.4	4.3	4.2	± 4.2	

LOT PX 403
Rev 2

CONTROL

 (Exp.) 2017-03-05
(YYYY - MM - DD)

		ABX Lysebio																							
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					L	CONTROL					N	CONTROL					H	TOLERANCES TOLERANCE					
		PENTRA						PENTRA						PENTRA											
		120	DX120	DX NEXUS				120	DX120	DX NEXUS				120	DX120	DX NEXUS					120	DX120	DX NEXUS		
		120 RETIC	DF120	DF NEXUS				120 RETIC	DF120	DF NEXUS				120 RETIC	DF120	DF NEXUS					120 RETIC	DF120	DF NEXUS		
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.5	2.5	2.5			± 0.4	7.7	7.7	7.7			± 1.0	19.3	19.3	19.3			± 2.2						
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.37	2.37	2.37			± 0.16	4.62	4.62	4.62			± 0.20	5.21	5.21	5.21			± 0.25						
	g/dl	6.9	6.9	6.9			± 0.4	13.7	13.7	13.7			± 0.5	16.1	16.1	16.1			± 0.6						
HB HGB	g/l	69	69	69			± 4	137	137	137			± 5	161	161	161			± 6						
	mmol/l	4.28	4.28	4.28			± 0.25	8.51	8.51	8.51			± 0.31	10.00	10.00	10.00			± 0.37						
HT HCT	%	19.7	19.7	19.7			± 1.5	37.9	37.9	37.9			± 2.0	45.3	45.3	45.3			± 2.5						
	l/l	0.197	0.197	0.197			± 0.015	0.379	0.379	0.379			± 0.020	0.453	0.453	0.453			± 0.025						
VGM MCV	µm ³ ·fl	83	83	83			± 5	82	82	82			± 5	87	87	87			± 5						
TGMH MCH	pg	29.1	29.1	29.1			± 2.0	29.7	29.7	29.7			± 2.0	30.9	30.9	30.9			± 2.5						
	fmol	1.81	1.81	1.81			± 0.12	1.84	1.84	1.84			± 0.12	1.92	1.92	1.92			± 0.16						
	g/dl	35.1	35.1	35.1			± 3.0	36.2	36.2	36.2			± 3.0	35.5	35.5	35.5			± 3.0						
CCMH MCHC	g/l	351	351	351			± 30	362	362	362			± 30	355	355	355			± 30						
	mmol/l	21.78	21.78	21.78			± 1.86	22.46	22.46	22.46			± 1.86	22.06	22.06	22.06			± 1.86						
IDR RDW	%	15.0	15.0	15.0			± 4.0	16.0	16.0	16.0			± 4.0	14.9	14.9	14.9			± 4.0						
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	72	72	72			± 20	244	244	244			± 30	485	485	485			± 50						
VPM MPV	µm ³ ·fl	10.2	10.2	10.2			± 2.0	10.0	10.0	10.0			± 2.0	10.0	10.0	10.0			± 2.0						
NEUT	#	1.49	1.53	1.53			± 0.35	4.24	4.37	4.37			± 0.90	13.70	13.90	13.90			± 1.90						
	%	59.4	61.3	61.3			± 10.0	55.1	56.7	56.7			± 10.0	71.2	72.2	72.2			± 10.0						
LYMPHO	#	0.68	0.62	0.62			± 0.33	2.49	2.33	2.33			± 0.70	2.68	2.47	2.47			± 1.50						
	%	27.1	24.7	24.7			± 12.0	32.3	30.2	30.2			± 8.0	13.9	12.8	12.8			± 8.0						
MONO	#	0.12	0.12	0.12			± 0.12	0.49	0.50	0.50			± 0.49	1.37	1.37	1.37			± 1.37						
	%	4.8	4.9	4.9			± 4.8	6.3	6.5	6.5			± 6.3	7.1	7.1	7.1			± 7.1						
EOS	#	0.16	0.17	0.17			± 0.16	0.25	0.28	0.28			± 0.25	1.02	1.04	1.04			± 1.02						
	%	6.2	6.6	6.6			± 6.2	3.3	3.6	3.6			± 3.3	5.3	5.4	5.4			± 5.3						
BASO	#	0.06	0.06	0.06			± 0.06	0.23	0.23	0.23			± 0.23	0.48	0.48	0.48			± 0.48						
	%	2.5	2.5	2.5			± 2.5	3.0	3.0	3.0			± 3.0	2.5	2.5	2.5			± 2.5						

1300027264-A
FRONT / RECTO
Ref: TEMP-0821 Rev.39