


LOT PX 400
Rev 2


CONTROL

 (Exp.) 2016-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio																		
			CONTROL					L	CONTROL					N	CONTROL					H	TOLERANCES TOLERANCE
			PENTRA						PENTRA						PENTRA						
			60 60C+ ES60	80 XL80	MS60	XLR	MSCRIP		60 60C+ ES60	80 XL80	MS60	XLR	MSCRIP		60 60C+ ES60	80 XL80	MS60	XLR	MSCRIP		
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.4	2.4	2.4	2.4	± 0.4	7.6	7.5	7.6	7.5	7.4	± 1.0	18.0	17.7	17.9	17.7	18.0	± 2.2	
GR	RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.41	2.37	2.36	2.37	2.36	± 0.16	4.62	4.64	4.58	4.64	4.60	± 0.20	5.06	5.12	5.06	5.12	5.08	± 0.25	
HB	HGB	g/dl	6.9	6.9	6.8	6.9	6.9	± 0.4	13.3	13.4	13.2	13.4	13.2	± 0.5	16.1	16.2	16.1	16.2	16.0	± 0.6	
		g/l	69	69	68	69	69	± 4	133	134	132	134	132	± 5	161	162	161	162	160	± 6	
		mmol/l	4.28	4.28	4.22	4.28	4.28	± 0.25	8.26	8.32	8.20	8.32	8.20	± 0.31	10.00	10.06	10.00	10.06	9.94	± 0.37	
HT	HCT	%	19.5	19.7	18.9	19.7	18.9	± 1.5	36.5	37.1	35.7	37.1	36.3	± 2.0	43.5	44.5	43.0	44.5	43.2	± 2.5	
		l/l	0.195	0.197	0.189	0.197	0.189	± 0.015	0.365	0.371	0.357	0.371	0.363	± 0.020	0.435	0.445	0.430	0.445	0.432	± 0.025	
VGM	MCV	µm ³ ; fl	81	83	80	83	80.0	± 5	79	80	78	80	79.0	± 5	86	87	85	87	85.0	± 5	
TGMH	MCH	pg	28.6	29.1	28.8	29.1	29.2	± 2.0	28.8	28.9	28.8	28.9	28.7	± 2.0	31.8	31.6	31.8	31.6	31.5	± 2.5	
		fmol	1.78	1.81	1.79	1.81	1.82	± 0.12	1.79	1.79	1.79	1.79	1.78	± 0.12	1.98	1.96	1.98	1.96	1.96	± 0.16	
CCMH	MCHC	g/dl	35.3	35.1	36.0	35.1	36.5	± 3.0	36.4	36.1	36.9	36.1	36.3	± 3.0	37.0	36.4	37.4	36.4	37.1	± 3.0	
		g/l	353	351	360	351	365	± 30	364	361	369	361	363	± 30	370	364	374	364	371	± 30	
		mmol/l	21.95	21.78	22.37	21.78	22.70	± 1.86	22.63	22.42	22.95	22.42	22.56	± 1.86	22.98	22.58	23.25	22.58	23.01	± 1.86	
IDR	RDW	%	13.3	13.0	12.2	13.0	11.7	± 4.0	13.0	13.5	12.2	13.5	12.5	± 4.0	12.8	13.0	11.8	13.0	12.0	± 4.0	
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	70	66	66	66	66	± 20	245	248	245	248	238	± 30	490	495	497	495	480	± 50	
VPM	MPV	µm ³ ; fl	9.9	10.4	10.1	10.4	9.4	± 2.0	9.5	9.9	9.8	9.9	9.2	± 2.0	8.7	9.1	9.0	9.1	8.3	± 2.0	
NEUT	#	#	1.37	1.41	1.41	1.41	1.38	± 0.35	4.07	4.05	4.10	4.05	3.97	± 0.90	12.60	12.39	12.62	12.39	12.69	± 1.90	
		%	57.0	58.6	58.7	58.6	57.5	± 10.0	53.5	54.0	54.0	54.0	53.7	± 10.0	70.0	70.0	70.5	70.0	70.5	± 10.0	
LYMPHO	#	#	0.70	0.67	0.67	0.67	0.72	± 0.33	2.51	2.45	2.49	2.45	2.48	± 0.70	2.81	2.81	2.81	2.81	2.84	± 1.50	
		%	29.0	28.0	28.0	28.0	29.8	± 12.0	33.0	32.7	32.8	32.7	33.5	± 8.0	15.6	15.9	15.7	15.9	15.8	± 8.0	
MONO	#	#	0.11	0.10	0.10	0.10	0.10	± 0.10	0.42	0.41	0.40	0.41	0.39	± 0.39	0.90	0.85	0.81	0.85	0.90	± 0.81	
		%	4.5	4.1	4.0	4.1	4.2	± 4.0	5.5	5.4	5.2	5.4	5.3	± 5.2	5.0	4.8	4.5	4.8	5.0	± 4.5	
EOS	#	#	0.14	0.14	0.14	0.14	0.12	± 0.12	0.34	0.34	0.34	0.34	0.30	± 0.30	0.90	0.89	0.90	0.89	0.81	± 0.81	
		%	6.0	6.0	5.8	6.0	5.0	± 5.0	4.5	4.5	4.5	4.5	4.0	± 4.0	5.0	5.0	5.0	5.0	4.5	± 4.5	
BASO	#	#	0.08	0.08	0.08	0.08	0.08	± 0.08	0.27	0.26	0.27	0.26	0.26	± 0.26	0.79	0.76	0.77	0.76	0.76	± 0.76	
		%	3.5	3.3	3.5	3.3	3.5	± 3.3	3.5	3.4	3.5	3.4	3.5	± 3.4	4.4	4.3	4.3	4.3	4.2	± 4.2	

LOT PX 400
Rev 2

CONTROL

 (Exp.) 2016-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio															
			CONTROL				L	CONTROL				N	CONTROL				H	TOLERANCES TOLERANCE
			PENTRA					PENTRA					PENTRA					
			120	DX120	DX NEXUS		120	DX120	DX NEXUS		120	DX120	DX NEXUS		120	DX120	DX NEXUS	
120 RETIC	DF120	DF NEXUS		120 RETIC	DF120	DF NEXUS		120 RETIC	DF120	DF NEXUS		120 RETIC	DF120	DF NEXUS				
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.4	2.4		± 0.4	7.8	7.8	7.8		± 1.0	18.6	18.6	18.6		± 2.2	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.40	2.40		± 0.16	4.63	4.63	4.63		± 0.20	5.14	5.14	5.14		± 0.25	
		g/dl	7.0	7.0	7.0		± 0.4	13.3	13.3	13.3		± 0.5	16.1	16.1	16.1		± 0.6	
HB	HGB	g/l	70	70	70		± 4	133	133	133		± 5	161	161	161		± 6	
		mmol/l	4.35	4.35	4.35		± 0.25	8.26	8.26	8.26		± 0.31	10.00	10.00	10.00		± 0.37	
HT	HCT	%	19.9	19.9	19.9		± 1.5	37.0	37.0	37.0		± 2.0	44.7	44.7	44.7		± 2.5	
		l/l	0.199	0.199	0.199		± 0.015	0.370	0.370	0.370		± 0.020	0.447	0.447	0.447		± 0.025	
VGM	MCV	µm ³ ; fl	83	83	83		± 5	80	80	80		± 5	87	87	87		± 5	
TGMH	MCH	pg	29.2	29.2	29.2		± 2.0	28.7	28.7	28.7		± 2.0	31.3	31.3	31.3		± 2.5	
		fmol	1.81	1.81	1.81		± 0.12	1.78	1.78	1.78		± 0.12	1.95	1.95	1.95		± 0.16	
		g/dl	35.1	35.1	35.1		± 3.0	35.9	35.9	35.9		± 3.0	36.0	36.0	36.0		± 3.0	
CCMH	MCHC	g/l	351	351	351		± 30	359	359	359		± 30	360	360	360		± 30	
		mmol/l	21.82	21.82	21.82		± 1.86	22.30	22.30	22.30		± 1.86	22.36	22.36	22.36		± 1.86	
IDR	RDW	%	15.7	15.7	15.7		± 4.0	17.0	17.0	17.0		± 4.0	14.9	14.9	14.9		± 4.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	70	70	70		± 20	250	250	250		± 30	490	490	490		± 50	
VPM	MPV	µm ³ ; fl	10.0	10.0	10.0		± 2.0	9.5	9.5	9.5		± 2.0	8.8	8.8	8.8		± 2.0	
		#	1.43	1.48	1.48		± 0.35	4.23	4.36	4.36		± 0.90	13.50	13.70	13.70		± 1.90	
		%	59.4	61.5	61.5		± 10.0	54.2	55.9	55.9		± 10.0	72.4	73.5	73.5		± 10.0	
		#	0.67	0.61	0.61		± 0.33	2.61	2.43	2.43		± 0.70	2.81	2.60	2.60		± 1.50	
		%	27.8	25.3	25.3		± 12.0	33.4	31.2	31.2		± 8.0	15.1	14.0	14.0		± 8.0	
		#	0.11	0.11	0.11		± 0.11	0.41	0.42	0.42		± 0.41	1.04	1.02	1.02		± 1.02	
		%	4.6	4.7	4.7		± 4.6	5.3	5.4	5.4		± 5.3	5.6	5.5	5.5		± 5.5	
		#	0.14	0.14	0.14		± 0.14	0.32	0.35	0.35		± 0.32	0.82	0.84	0.84		± 0.82	
		%	5.7	6.0	6.0		± 5.7	4.1	4.5	4.5		± 4.1	4.4	4.5	4.5		± 4.4	
		#	0.06	0.06	0.06		± 0.06	0.23	0.23	0.23		± 0.23	0.47	0.47	0.47		± 0.47	
		%	2.5	2.5	2.5		± 2.5	3.0	3.0	3.0		± 3.0	2.5	2.5	2.5		± 2.5	

Ref: TEMP-0821 Rev.37 FRONT / RECTO 1300019740-A