



LOT PX 113

CONTROL

(Exp.) **2014-01-05**
 (YYYY-MM-DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse																				
				CONTROL				L			CONTROL				N			CONTROL				H		
				PENTRA							PENTRA							PENTRA						
				60 60C+ ES60	80 XL80			TOLERANCES TOLERANCE	60 60C+ ES60	80 XL80			TOLERANCES TOLERANCE	60 60C+ ES60	80 XL80			TOLERANCES TOLERANCE	60 60C+ ES60	80 XL80			TOLERANCES TOLERANCE	
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.5				± 0.4	7.4	7.3				± 1.0	17.7	17.5				± 2.2				
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.39	2.33				± 0.12	4.58	4.56				± 0.15	5.15	5.15				± 0.20				
HB	HGB	g/dl	6.8	6.7				± 0.4	13.4	13.4				± 0.5	16.3	16.2				± 0.6				
		g/l	68	67				± 4	134	134				± 5	163	162				± 6				
		mmol/l	4.22	4.16				± 0.25	8.32	8.32				± 0.31	10.12	10.06				± 0.37				
HT	HCT	%	18.9	19.1				± 1.5	37.1	37.4				± 2.0	44.8	45.3				± 2.5				
		l/l	0.189	0.191				± 0.015	0.371	0.374				± 0.020	0.448	0.453				± 0.025				
VGM	MCV	µm ³ ; fl	79	82				± 5	81	82				± 5	87	88				± 5				
TGMH	MCH	pg	28.5	28.8				± 2.0	29.3	29.4				± 2.0	31.7	31.5				± 2.5				
		fmol	1.77	1.79				± 0.12	1.82	1.82				± 0.12	1.97	1.95				± 0.16				
CCMH	MCHC	g/dl	36.0	35.1				± 3.0	36.1	35.8				± 3.0	36.4	35.7				± 3.0				
		g/l	360	351				± 30	361	358				± 30	364	357				± 30				
		mmol/l	22.37	21.78				± 1.86	22.43	22.25				± 1.86	22.59	22.20				± 1.86				
IDR	RDW	%	13.0	14.0				± 4.0	12.5	13.8				± 4.0	11.6	13.3				± 4.0				
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	72	70				± 20	258	260				± 30	505	515				± 50				
VPM	MPV	µm ³ ; fl	9.0	9.2				± 2.0	8.7	9.0				± 2.0	8.8	8.9				± 2.0				
NEUT	#	1.36	1.39				± 0.30	3.96	3.83				± 0.80	12.74	12.43				± 1.90					
	%	56.5	55.5				± 10.0	53.5	52.5				± 10.0	72.0	71.0				± 10.0					
LYMPHO	#	0.73	0.79				± 0.25	2.59	2.63				± 0.60	2.74	2.80				± 1.50					
	%	30.5	31.5				± 8.0	35.0	36.0				± 8.0	15.5	16.0				± 8.0					
MONO	#	0.08	0.08				± 0.08	0.33	0.30				± 0.30	0.72	0.72				± 0.72					
	%	3.0	3.0				± 3.0	4.5	4.0				± 4.0	4.0	4.0				± 4.0					
EOS	#	0.16	0.16				± 0.11	0.26	0.29				± 0.19	0.71	0.79				± 0.45					
	%	6.5	6.5				± 4.5	3.5	4.0				± 2.5	4.0	4.5				± 2.5					
BASO	#	0.08	0.09				± 0.07	0.26	0.26				± 0.23	0.80	0.79				± 0.45					
	%	3.5	3.5				± 3.0	3.5	3.5				± 3.0	4.5	4.5				± 2.5					



ABX Difftrol

LOT PX 113

CONTROL



(Exp.) 2014-01-05
 (YYYY-MM-DD)

Rev 2

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse																		
		CONTROL				L	TOLERANCES TOLERANCE	CONTROL				N	TOLERANCES TOLERANCE	CONTROL				H	TOLERANCES TOLERANCE	
		PENTRA				PENTRA				PENTRA				PENTRA						
		120 120 RETIC	DX120 DF120			120 120 RETIC		DX120 DF120			120 120 RETIC	DX120 DF120				120 120 RETIC	DX120 DF120			
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l					± 0.4						± 1.0					± 2.2			
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l					± 0.12						± 0.15					± 0.20			
HB HGB	g/dl					± 0.4						± 0.5					± 0.6			
	g/l					± 4						± 5					± 6			
	mmol/l					± 0.25						± 0.31					± 0.37			
HT HCT	%					± 1.5						± 2.0					± 2.5			
	l/l					± 0.015						± 0.020					± 0.025			
VGM MCV	µm ³ ·fl					± 5						± 5					± 5			
TGMH MCH	pg					± 2.0						± 2.0					± 2.5			
	fmol					± 0.12						± 0.12					± 0.16			
CCMH MCHC	g/dl					± 3.0						± 3.0					± 3.0			
	g/l					± 30						± 30					± 30			
	mmol/l					± 1.86						± 1.86					± 1.86			
IDR RDW	%					± 4.0						± 4.0					± 4.0			
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l					± 20						± 30					± 50			
VPM MPV	µm ³ ·fl					± 2.0						± 2.0					± 2.0			
NEUT	#					± 0.30						± 0.80					± 1.90			
	%					± 10.0						± 10.0					± 10.0			
LYMPHO	#					± 0.25						± 0.60					± 1.50			
	%					± 8.0						± 8.0					± 8.0			
MONO	#					± 0.08						± 0.30					± 0.72			
	%					± 3.0						± 4.0					± 4.0			
EOS	#					± 0.11						± 0.19					± 0.45			
	%					± 4.5						± 2.5					± 2.5			
BASO	#					± 0.07						± 0.23					± 0.45			
	%					± 3.0						± 3.0					± 2.5			

Ref: TEMP-0821 Rev.33 FRONT / RECTO 9930080-B