


LOT
PX 053
CONTROL

(Exp.) 2013-07-05

(YYYY-MM-DD)

Rev 3

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse																	
		CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE
		PENTRA						PENTRA						PENTRA					
		60 60C+ ES60	80 XL80					60 60C+ ES60	80 XL80					60 60C+ ES60	80 XL80				
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.5				± 0.4	7.4	7.4				± 1.0	17.2	17.2				± 2.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.48	2.44				± 0.12	4.56	4.57				± 0.15	5.02	5.05				± 0.20
HB HGB	g/dl	7.0	7.0				± 0.4	13.7	13.8				± 0.5	15.9	15.9				± 0.6
	g/l	70	70				± 4	137	138				± 5	159	159				± 6
	mmol/l	4.35	4.35				± 0.25	8.51	8.57				± 0.31	9.87	9.87				± 0.37
HT HCT	%	19.8	20.3				± 1.5	37.8	38.4				± 2.0	43.7	43.9				± 2.5
	l/l	0.198	0.203				± 0.015	0.378	0.384				± 0.020	0.437	0.439				± 0.025
VGM MCV	µm ³ ; fl	80	83				± 5	83	84				± 5	87	87				± 5
TGMH MCH	pg	28.2	28.7				± 2.0	30.0	30.2				± 2.0	31.7	31.5				± 2.5
	fmol	1.75	1.78				± 0.12	1.87	1.88				± 0.12	1.97	1.96				± 0.16
CCMH MCHC	g/dl	35.3	34.6				± 3.0	36.2	35.9				± 3.0	36.4	36.2				± 3.0
	g/l	353	346				± 30	362	359				± 30	364	362				± 30
	mmol/l	21.91	21.46				± 1.86	22.48	22.32				± 1.86	22.61	22.47				± 1.86
IDR RDW	%	13.2	13.6				± 4.0	12.5	13.6				± 4.0	12.3	13.2				± 4.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	72	70				± 20	258	260				± 30	505	515				± 50
VPM MPV	µm ³ ; fl	9.2	9.3				± 2.0	8.9	9.1				± 2.0	8.8	9.0				± 2.0
NEUT	#	1.49	1.60				± 0.30	4.37	4.40				± 0.80	11.87	11.78				± 1.90
	%	62.0	64.0				± 10.0	59.0	59.5				± 10.0	69.0	68.5				± 10.0
LYMPHO	#	0.56	0.55				± 0.30	2.18	2.11				± 0.70	2.92	2.92				± 1.50
	%	23.5	22.0				± 9.0	29.5	28.5				± 8.0	17.0	17.0				± 8.0
MONO	#	0.13	0.11				± 0.11	0.33	0.33				± 0.33	0.77	0.77				± 0.77
	%	5.5	4.5				± 4.5	4.5	4.5				± 4.5	4.5	4.5				± 4.5
EOS	#	0.13	0.15				± 0.13	0.27	0.30				± 0.27	0.86	0.95				± 0.82
	%	5.5	6.0				± 5.0	3.5	4.0				± 3.5	5.0	5.5				± 4.5
BASO	#	0.08	0.09				± 0.05	0.26	0.26				± 0.16	0.77	0.77				± 0.36
	%	3.5	3.5				± 2.0	3.5	3.5				± 2.0	4.5	4.5				± 2.0



LOT PX 053

CONTROL



(Exp.) 2013-07-05
 (YYYY - MM - DD)

Rev 3

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse																		
		CONTROL				L	TOLERANCES TOLERANCE	CONTROL				N	TOLERANCES TOLERANCE	CONTROL				H	TOLERANCES TOLERANCE	
		PENTRA						PENTRA						PENTRA						
		120 120 RETIC	DX120 DF120						120 120 RETIC	DX120 DF120						120 120 RETIC	DX120 DF120			
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.5	2.5				± 0.4	7.8	7.8				± 1.0	18.2	18.2				± 2.2	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.48	2.48				± 0.12	4.60	4.60				± 0.15	5.14	5.14				± 0.20	
HB HGB	g/dl	7.1	7.1				± 0.4	13.8	13.8				± 0.5	15.9	15.9				± 0.6	
	g/l	71	71				± 4	138	138				± 5	159	159				± 6	
	mmol/l	4.41	4.41				± 0.25	8.57	8.57				± 0.31	9.87	9.87				± 0.37	
HT HCT	%	20.3	20.3				± 1.5	38.6	38.6				± 2.0	44.7	44.7				± 2.5	
	l/l	0.203	0.203				± 0.015	0.386	0.386				± 0.020	0.447	0.447				± 0.025	
VGM MCV	µm ³ ; fl	82	82				± 5	84	84				± 5	87	87				± 5	
TGMH MCH	pg	28.6	28.6				± 2.0	30.0	30.0				± 2.0	30.9	30.9				± 2.5	
	fmol	1.78	1.78				± 0.12	1.86	1.86				± 0.12	1.92	1.92				± 0.16	
CCMH MCHC	g/dl	34.9	34.9				± 3.0	35.7	35.7				± 3.0	35.6	35.6				± 3.0	
	g/l	349	349				± 30	357	357				± 30	356	356				± 30	
	mmol/l	21.68	21.68				± 1.86	22.18	22.18				± 1.86	22.08	22.08				± 1.86	
IDR RDW	%	16.0	16.0				± 4.0	16.0	16.0				± 4.0	15.2	15.2				± 4.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	74	74				± 20	258	258				± 30	508	508				± 50	
VPM MPV	µm ³ ; fl	9.1	9.1				± 2.0	9.0	9.0				± 2.0	8.9	8.9				± 2.0	
NEUT	#	1.60	1.55				± 0.30	4.60	4.60				± 0.80	12.69	12.65				± 1.90	
	%	64.0	62.0				± 10.0	59.0	59.0				± 10.0	69.7	69.5				± 10.0	
LYMPHO	#	0.60	0.56				± 0.30	2.34	2.18				± 0.70	3.24	3.00				± 1.50	
	%	24.0	22.5				± 9.0	30.0	28.0				± 8.0	17.8	16.5				± 8.0	
MONO	#	0.11	0.19				± 0.11	0.39	0.51				± 0.33	1.00	1.18				± 0.77	
	%	4.5	7.5				± 4.5	5.0	6.5				± 4.5	5.5	6.5				± 4.5	
EOS	#	0.13	0.13				± 0.13	0.27	0.27				± 0.27	0.82	0.82				± 0.82	
	%	5.0	5.0				± 5.0	3.5	3.5				± 3.5	4.5	4.5				± 4.5	
BASO	#	0.06	0.08				± 0.05	0.20	0.23				± 0.16	0.46	0.55				± 0.36	
	%	2.5	3.0				± 2.0	2.5	3.0				± 2.0	2.5	3.0				± 2.0	

Ref: TEMP-0821 Rev.32 FRONT / RECTO 9930080-A