



LOT PX 073

CONTROL

(Exp.) **2013-09-05**
 (YYYY-MM-DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse																				
				CONTROL				L	TOLERANCES TOLERANCE		CONTROL				N	TOLERANCES TOLERANCE		CONTROL				H	TOLERANCES TOLERANCE	
				PENTRA				PENTRA				PENTRA												
				60 60C+ ES60	80 XL80			60 60C+ ES60	80 XL80			60 60C+ ES60	80 XL80			60 60C+ ES60	80 XL80							
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.4	2.4				± 0.4	7.3	7.3				± 1.0	17.8	17.5				± 2.2				
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.42	2.37				± 0.12	4.58	4.58				± 0.15	5.10	5.11				± 0.20				
HB	HGB	g/dl	6.8	6.9				± 0.4	13.2	13.3				± 0.5	16.0	16.0				± 0.6				
		g/l	68	69				± 4	132	133				± 5	160	160				± 6				
		mmol/l	4.22	4.28				± 0.25	8.20	8.26				± 0.31	9.94	9.94				± 0.37				
HT	HCT	%	19.1	19.7				± 1.5	36.2	36.6				± 2.0	44.4	44.5				± 2.5				
		l/l	0.191	0.197				± 0.015	0.362	0.366				± 0.020	0.444	0.445				± 0.025				
VGM	MCV	µm ³ ; fl	79	83				± 5	79	80				± 5	87	87				± 5				
TGMH	MCH	pg	28.1	29.1				± 2.0	28.8	29.0				± 2.0	31.4	31.3				± 2.5				
		fmol	1.74	1.81				± 0.12	1.79	1.80				± 0.12	1.95	1.94				± 0.16				
CCMH	MCHC	g/dl	35.6	35.1				± 3.0	36.5	36.3				± 3.0	36.1	36.0				± 3.0				
		g/l	356	351				± 30	365	363				± 30	361	360				± 30				
		mmol/l	22.09	21.78				± 1.86	22.66	22.54				± 1.86	22.39	22.35				± 1.86				
IDR	RDW	%	12.7	13.7				± 4.0	12.5	13.6				± 4.0	12.3	13.4				± 4.0				
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	75	72				± 20	250	245				± 30	510	505				± 50				
VPM	MPV	µm ³ ; fl	9.2	9.4				± 2.0	9.4	9.6				± 2.0	9.3	9.6				± 2.0				
NEUT	#		1.36	1.36				± 0.30	4.23	4.25				± 0.80	12.82	12.60				± 1.90				
	%		56.7	56.7				± 10.0	57.9	58.2				± 10.0	72.0	72.0				± 10.0				
LYMPHO	#		0.73	0.72				± 0.25	2.28	2.21				± 0.65	2.58	2.54				± 1.50				
	%		30.5	30.0				± 8.0	31.3	30.3				± 8.0	14.5	14.5				± 8.0				
MONO	#		0.08	0.08				± 0.08	0.30	0.30				± 0.30	0.71	0.62				± 0.62				
	%		3.5	3.3				± 3.3	4.0	4.0				± 4.0	4.0	3.5				± 3.5				
EOS	#		0.14	0.16				± 0.14	0.27	0.29				± 0.27	0.89	0.96				± 0.83				
	%		6.0	6.5				± 5.5	3.5	4.0				± 3.5	5.0	5.5				± 4.5				
BASO	#		0.08	0.08				± 0.08	0.24	0.26				± 0.23	0.80	0.79				± 0.46				
	%		3.3	3.5				± 3.0	3.3	3.5				± 3.0	4.5	4.5				± 2.5				



ABX Difftrol

LOT PX 073

CONTROL



(Exp.) 2013-09-05
 (YYYY-MM-DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lyse																	
			CONTROL				L	TOLERANCES TOLERANCE	CONTROL				N	TOLERANCES TOLERANCE	CONTROL				H	TOLERANCES TOLERANCE
			PENTRA						PENTRA						PENTRA					
			120 120 RETIC	DX120 DF120			120 120 RETIC		DX120 DF120			120 120 RETIC	DX120 DF120				120 120 RETIC	DX120 DF120		
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l					± 0.4						± 1.0						± 2.2	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l					± 0.12						± 0.15						± 0.20	
HB	HGB	g/dl					± 0.4						± 0.5						± 0.6	
		g/l					± 4					± 5					± 6			
		mmol/l					± 0.25					± 0.31					± 0.37			
HT	HCT	%					± 1.5					± 2.0						± 2.5		
		l/l					± 0.015					± 0.020					± 0.025			
VGM	MCV	µm ³ ·fl					± 5					± 5						± 5		
TGMH	MCH	pg					± 2.0					± 2.0						± 2.5		
		fmol					± 0.12					± 0.12					± 0.16			
CCMH	MCHC	g/dl					± 3.0					± 3.0						± 3.0		
		g/l					± 30					± 30					± 30			
		mmol/l					± 1.86					± 1.86					± 1.86			
IDR	RDW	%					± 4.0					± 4.0					± 4.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l					± 20					± 30						± 50		
VPM	MPV	µm ³ ·fl					± 2.0					± 2.0						± 2.0		
NEUT	#						± 0.35					± 0.90						± 1.90		
	%						± 10.0					± 10.0						± 10.0		
LYMPHO	#						± 0.25					± 0.70						± 1.50		
	%						± 8.0					± 8.0						± 8.0		
MONO	#						± 0.03					± 0.07						± 0.18		
	%						± 1.0					± 1.0						± 1.0		
EOS	#						± 0.15					± 0.24						± 0.70		
	%						± 4.0					± 3.0						± 3.0		
BASO	#						± 0.06					± 0.25						± 0.60		
	%						± 2.5					± 3.0						± 3.0		

Ref: TEMP-0821 Rev.33 FRONT / RECTO 9930080-B