


**LOT** PX 427  
Rev 2


**CONTROL**

 (Exp.) 2021-03-05  
(YYYY - MM - DD)

		ABX Lysebio																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				L	TOLERANCES TOLERANCE	CONTROL				N	TOLERANCES TOLERANCE	CONTROL				H	TOLERANCES TOLERANCE
		PENTRA						PENTRA						PENTRA					
		80 XL80	XLR					80 XL80	XLR					80 XL80	XLR				
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	2.7	2.7			± 0.4	7.8	7.8			± 1.0	17.8	17.8			± 2.2			
GR RBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	2.29	2.29			± 0.16	4.50	4.50			± 0.20	5.11	5.11			± 0.25			
HB HGB	g/dl	6.6	6.6			± 0.4	13.7	13.7			± 0.5	16.8	16.8			± 0.6			
	g/l	66	66			± 4	137	137			± 5	168	168			± 6			
	mmol/l	4.10	4.10			± 0.25	8.51	8.51			± 0.31	10.43	10.43			± 0.37			
HT HCT	%	19.2	19.2			± 1.5	37.8	37.8			± 2.0	46.5	46.5			± 2.5			
	l/l	0.192	0.192			± 0.015	0.378	0.378			± 0.020	0.465	0.465			± 0.025			
VGM MCV	µm <sup>3</sup> fl	84	84			± 5	84	84			± 5	91	91			± 5			
TGMH MCH	pg	28.8	28.8			± 2.0	30.4	30.4			± 2.0	32.9	32.9			± 2.5			
	fmol	1.79	1.79			± 0.12	1.89	1.89			± 0.12	2.04	2.04			± 0.16			
CCMH MCHC	g/dl	34.3	34.3			± 3.0	36.2	36.2			± 3.0	36.1	36.1			± 3.0			
	g/l	343	343			± 30	362	362			± 30	361	361			± 30			
	mmol/l	21.31	21.31			± 1.86	22.51	22.51			± 1.86	22.44	22.44			± 1.86			
IDR RDW	%	13.7	13.3			± 4.0	14.0	13.1			± 4.0	13.0	12.5			± 4.0			
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	68	68			± 20	252	252			± 30	490	490			± 50			
VPM MPV	µm <sup>3</sup> fl	9.3	9.3			± 2.0	9.1	9.1			± 2.0	9.5	9.5			± 2.0			
NEUT	#	1.41	1.41			± 0.35	4.02	4.02			± 0.90	11.71	11.71			± 1.90			
	%	52.3	52.3			± 10.0	51.5	51.5			± 10.0	65.8	65.8			± 10.0			
LYMPHO	#	0.77	0.77			± 0.33	2.36	2.36			± 0.70	2.74	2.74			± 1.50			
	%	28.7	28.7			± 12.0	30.2	30.2			± 8.0	15.4	15.4			± 8.0			
MONO	#	0.23	0.23			± 0.23	0.76	0.76			± 0.76	1.66	1.66			± 1.66			
	%	8.4	8.4			± 8.4	9.7	9.7			± 9.7	9.3	9.3			± 9.3			
EOS	#	0.19	0.19			± 0.19	0.37	0.37			± 0.37	0.87	0.87			± 0.87			
	%	7.2	7.2			± 7.2	4.8	4.8			± 4.8	4.9	4.9			± 4.9			
BASO	#	0.09	0.09			± 0.09	0.30	0.30			± 0.30	0.82	0.82			± 0.82			
	%	3.4	3.4			± 3.4	3.8	3.8			± 3.8	4.6	4.6			± 4.6			

**LOT** PX 427  
Rev 2

**CONTROL**

 (Exp.) 2021-03-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio															TOLERANCES TOLERANCE
			CONTROL			L	CONTROL			N	CONTROL			H	TOLERANCES TOLERANCE			
			PENTRA			TOLERANCES TOLERANCE	PENTRA			TOLERANCES TOLERANCE	PENTRA			TOLERANCES TOLERANCE				
			60 60C+ ES60	MS60	MSCRP		60 60C+ ES60	MS60	MSCRP		60 60C+ ES60	MS60	MSCRP					
GB	WBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	2.6	2.6	2.7		± 0.4	7.7	7.8	7.8		± 1.0	18.5	18.2	18.4		± 2.2	
GR	RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	2.34	2.31	2.33		± 0.16	4.51	4.50	4.51		± 0.20	5.08	5.04	5.14		± 0.25	
		g/dl	6.6	6.5	6.6		± 0.4	13.6	13.4	13.5		± 0.5	16.7	16.5	16.6		± 0.6	
HB	HGB	g/l	66	65	66		± 4	136	134	135		± 5	167	165	166		± 6	
		mmol/l	4.10	4.04	4.10		± 0.25	8.45	8.32	8.38		± 0.31	10.37	10.25	10.31		± 0.37	
HT	HCT	%	18.7	18.5	18.5		± 1.5	37.4	36.9	37.0		± 2.0	45.7	44.9	45.2		± 2.5	
		l/l	0.187	0.185	0.185		± 0.015	0.374	0.369	0.370		± 0.020	0.457	0.449	0.452		± 0.025	
VGM	MCV	µm <sup>3</sup> ; fl	80	80	79.5		± 5	83	82	82.0		± 5	90	89	88.0		± 5	
TGMH	MCH	pg	28.2	28.1	28.3		± 2.0	30.2	29.8	29.9		± 2.0	32.9	32.7	32.3		± 2.5	
		fmol	1.75	1.75	1.76		± 0.12	1.87	1.85	1.86		± 0.12	2.04	2.03	2.01		± 0.16	
		g/dl	35.3	35.2	35.6		± 3.0	36.3	36.3	36.5		± 3.0	36.5	36.8	36.7		± 3.0	
CCMH	MCHC	g/l	353	352	356		± 30	363	363	365		± 30	365	368	367		± 30	
		mmol/l	21.89	21.84	22.13		± 1.86	22.56	22.55	22.67		± 1.86	22.68	22.84	22.79		± 1.86	
IDR	RDW	%	13.5	12.5	13.0		± 4.0	12.5	12.5	12.8		± 4.0	12.0	12.0	12.2		± 4.0	
PLAQ.	PLTS	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	72	70	71		± 20	250	251	252		± 30	491	490	489		± 50	
VPM	MPV	µm <sup>3</sup> ; fl	8.5	8.6	8.4		± 2.0	8.5	8.3	8.1		± 2.0	9.0	8.8	8.5		± 2.0	
		#	1.46	1.42	1.51		± 0.35	4.27	4.20	4.33		± 0.90	12.95	12.38	13.06		± 1.90	
		%	56.2	54.6	56.0		± 10.0	55.5	53.9	55.5		± 10.0	70.0	68.0	71.0		± 10.0	
		#	0.77	0.79	0.78		± 0.33	2.51	2.50	2.46		± 0.70	2.96	2.84	2.76		± 1.50	
		%	29.8	30.5	29.0		± 12.0	32.6	32.1	31.5		± 8.0	16.0	15.6	15.0		± 8.0	
		#	0.11	0.14	0.14		± 0.11	0.32	0.51	0.39		± 0.32	0.81	1.22	0.92		± 0.81	
		%	4.2	5.3	5.0		± 4.2	4.2	6.6	5.0		± 4.2	4.4	6.7	5.0		± 4.4	
		#	0.16	0.16	0.16		± 0.16	0.29	0.29	0.31		± 0.29	0.96	0.95	0.92		± 0.92	
		%	6.3	6.1	6.0		± 6.0	3.8	3.7	4.0		± 3.7	5.2	5.2	5.0		± 5.0	
		#	0.09	0.09	0.11		± 0.09	0.30	0.29	0.31		± 0.29	0.81	0.82	0.74		± 0.74	
		%	3.5	3.5	4.0		± 3.5	3.9	3.7	4.0		± 3.7	4.4	4.5	4.0		± 4.0	

1300040596 FRONT / RECTO Ref: TEMP-0821 Rev.44