



Yumizen G PT Reco 10 for Yumizen G Line



D 4 9 2 0 4 2 1 3 6 3 7 6 0 4 2 4 D



D 5 1 0 0 1 0 5 5 0 0 1 5 0 3 3 0 D



D 6 1 9 6 2 5 0 2 4 0 1 0 5 1 0 8 D



D 7 1 4 3 D

Reference No.:		1300036376	
Lot No.:		920421 (rev. 1)	
Expiry date:		2024-04-30	
Coagulometers		ISI	MASTER CURVE
Optical	Yumizen G Line		%
		100	10,5
		50	15,0
		33	19,6
		25	24,0
		MNPT	10,5
Important notice			
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p> <p>For Yumizen G200, Yumizen G400 and Yumizen G400 DDi, the Min. Time must be set to 7 seconds and the Lag time to 5 seconds in the PT setup. Please refer to the Yumizen G PT Reco 10 reagent notice.</p>			