

# Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



# Yumizen G FIB 5 for Yumizen G400 DDi/G200



<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>920408 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-04-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>g/L</b>	<b>sec</b>
<b>Optical</b>	Yumizen G1500/G1550/G800	<b>4</b>	<b>7,6</b>
		<b>3</b>	<b>9,3</b>
		<b>2</b>	<b>12,3</b>
		<b>1</b>	<b>19,8</b>
		<b>Dilution</b>	<b>10x</b>

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>920408 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-04-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>g/L</b>	<b>sec</b>
<b>Optical</b>	Yumizen G400 DDi/G200	<b>4</b>	<b>8,6</b>
		<b>3</b>	<b>10,3</b>
		<b>2</b>	<b>13,8</b>
		<b>1</b>	<b>23,2</b>
		<b>Dilution</b>	<b>10x</b>