

# Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4920317363910323D



D 5400264030020002D



D 6001230100054005D



D 7002600250135142D

Reference No.:		1300036391	
Lot No.:		920317 (rev. 1)	
Expiry date:		2023-03-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0,2640
		3,00	0,2000
		2,00	0,1230
		1,00	0,0540
		0,50	0,0260
		0,25	0,0135
Slope value		0,06	

# Yumizen G DDi 2 for YUMIZEN G800



D 4920317363910323D



D 5400246530018852D



D 6001150100049505D



D 7002350250115177D

Reference No.:		1300036391	
Lot No.:		920317 (rev. 1)	
Expiry date:		2023-03-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0,2465
		3,00	0,1885
		2,00	0,1150
		1,00	0,0495
		0,50	0,0235
		0,25	0,0115
Slope value		0,06	



# Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

<b>Reference No.:</b>		<b>1300036391</b>	
<b>Lot No.:</b>		<b>920317 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2023-03-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>µg(FEU)/mL</b>	<b>ΔOD</b>
<b>Optical</b>	<b>Yumizen G400 DDi/G200</b>	<b>4,00</b>	<b>0,1009</b>
		<b>3,00</b>	<b>0,0785</b>
		<b>2,00</b>	<b>0,0469</b>
		<b>1,00</b>	<b>0,0192</b>
		<b>0,50</b>	<b>0,0090</b>
		<b>0,25</b>	<b>0,0045</b>
		<b>Slope value</b>	<b>0,02</b>
<b>Important notice</b>			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type <b>Lin-Lin p-p</b> is selected in <b>Calibration &gt; D-Dimer &gt; ug FEU/ml</b>.          For more information, please contact your local HORIBA Medical representative.</p>			