



Yumizen G PT 5 for Yumizen G Line



D 4 9 0 0 8 1 1 3 6 3 3 8 0 8 2 2 D



D 5 1 0 0 1 3 1 5 0 0 1 7 1 3 3 0 D



D 6 2 1 4 2 5 0 2 5 3 1 3 1 1 2 7 D



D 7 1 4 1 D

Reference No.:		1300036338		
Lot No.:		900811 (rev. 3)		
Expiry date:		2022-08-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,27	%
		100		13,1
		50		17,1
		33		21,4
		25		25,3
			MNPT	13,1
Important notice				
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be improved by performing it at 37°C. 				