

Yumizen G PT 5 for Yumizen G Line



Reference No.:		1300036338			
Lot No.:		911003 (rev. 2)			
Expiry date:		2023-10-31			
Coagulometers		ISI	MASTER CURVE		
Optical	Yumizen G Line		%	sec	
		1,18	100	12,3	
			50	16,3	
			33	20,5	
			25	24,5	
		MNPT	12,3		
Important notice					
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be improved by performing it at 37°C. 					