



Yumizen G PT Reco 10 for Yumizen G Line



D 4910620363760623D



D 5100101500131330D



D 6162250191101105D



D 7130D

Reference No.:		1300036376	
Lot No.:		910620 (rev. 2)	
Expiry date:		2023-06-30	
Coagulometers		ISI	MASTER CURVE
Optical	Yumizen G Line		%
		100	10,1
		50	13,1
		33	16,2
		25	19,1
		MNPT	10,1
Important notice			
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p> <p>For Yumizen G200, Yumizen G400 and Yumizen G400 DDi, the Min. Time must be set to 7 seconds and the Lag time to 5 seconds in the PT setup. Please refer to the Yumizen G PT Reco 10 reagent notice.</p>			