



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4 9 1 0 3 1 8 4 6 4 1 2 0 3 2 3 D



D 5 1 1 2 1 5 2 0 6 6 1 1 1 0 8 6 D



D 6 1 2 9 1 1 6 1 5 6 0 5 5 0 9 2 D



D 7 0 9 3 1 3 9 0 8 8 1 2 0 0 7 3 D



D 8 1 2 1 0 8 1 1 2 1 3 3 9 4 5 9 D



D 9 3 0 8 4 1 6 2 0 0 3 0 0 2 0 1 D



D 0 3 0 1 3 1 8 D

Reference No.:		1300036412	
Lot No:		910318 (rev. 4)	
Expiry date:		2023-03-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	11,2-15,2
	Prothrombin %	%	66-111
	Prothrombin INR	INR	0,86-1,29
PT Liq (*)	Prothrombin Time	sec	11,6-15,6 (*)
	Prothrombin %	%	55-92 (*)
	Prothrombin INR	INR	0,93-1,39 (*)
PT R	Prothrombin Time	sec	8,8-12,0
	Prothrombin %	%	73-121
	Prothrombin INR	INR	0,81-1,21
APTT	Partial Thromboplastin Time	sec	33,9-45,9
APTT Liq	Partial Thromboplastin Time	sec	30,8-41,6
FIB Cl	Fibrinogen	g/L	2,0-3,0
TT	Thrombin Time	sec	20,1-30,1
AT	Antithrombin	%	86-116
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			

(*) Warning! This range has to be used only with Yumizen G PT Liq 4 lots : 920315 and higher except for the Yumizen G PT Liq 4 lots 920611 & 920705. For the other lots, including 920611 & 920705, please refer to the previous revision of values sheet.



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 1 0 3 1 8 5 6 4 1 2 0 3 2 3 D



D 5 1 4 9 2 2 4 0 3 5 0 5 8 1 3 1 D



D 6 1 9 6 1 5 4 2 3 1 0 2 8 0 4 7 D



D 7 1 4 1 2 1 2 1 1 1 6 6 0 3 9 D



D 8 0 6 5 1 0 9 1 6 3 4 7 8 6 4 7 D



D 9 4 4 4 6 0 0 0 8 0 1 8 0 0 0 0 D



D 0 9 9 9 3 5 6 D

Reference No.:		1300036412	
Lot No:		910318 (rev. 4)	
Expiry date:		2023-03-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	14,9-22,4
	Prothrombin %	%	35-58
	Prothrombin INR	INR	1,31-1,96
PT Liq (*)	Prothrombin Time	sec	15,4-23,1 (*)
	Prothrombin %	%	28-47 (*)
	Prothrombin INR	INR	1,41-2,12 (*)
PT R	Prothrombin Time	sec	11,1-16,6
	Prothrombin %	%	39-65
	Prothrombin INR	INR	1,09-1,63
APTT	Partial Thromboplastin Time	sec	47,8-64,7
APTT Liq	Partial Thromboplastin Time	sec	44,4-60,0
FIB Cl	Fibrinogen	g/L	0,8-1,8
TT	Thrombin Time	sec	NA
AT	Antithrombin	%	41-68

Important notice

The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.

The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.

(*) Warning! This range has to be used only with Yumizen G PT Liq 4 lots : 920315 and higher except for the Yumizen G PT Liq 4 lots 920611 & 920705. For the other lots, including 920611 & 920705, please refer to the previous revision of values sheet.