



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4 9 1 0 4 3 4 4 6 4 1 2 0 4 2 3 D



D 5 1 1 2 1 5 2 0 6 6 1 1 0 0 8 6 D



D 6 1 2 9 1 2 6 1 7 1 0 6 1 1 0 1 D



D 7 0 9 2 1 3 8 0 8 6 1 1 6 0 7 9 D



D 8 1 3 1 0 7 8 1 1 8 3 2 2 4 3 6 D



D 9 3 0 0 4 0 6 2 1 0 3 1 0 2 0 1 D



D 0 3 0 1 3 0 3 D

Reference No.:		1300036412	
Lot No:		910434 (rev. 3)	
Expiry date:		2023-04-30	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	11,2-15,2
	Prothrombin %	%	66-110
	Prothrombin INR	INR	0,86-1,29
PT Liq (*)	Prothrombin Time	sec	12,6-17,1 (*)
	Prothrombin %	%	61-101 (*)
	Prothrombin INR	INR	0,92-1,38 (*)
PT R	Prothrombin Time	sec	8,6-11,6
	Prothrombin %	%	79-131
	Prothrombin INR	INR	0,78-1,18
APTT	Partial Thromboplastin Time	sec	32,2-43,6
APTT Liq	Partial Thromboplastin Time	sec	30,0-40,6
FIB Cl	Fibrinogen	g/L	2,1-3,1
TT	Thrombin Time	sec	20,1-30,1
AT	Antithrombin	%	85-114
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			

(*) Warning! This range has to be used only with Yumizen G PT Liq 4 lots : 920315 and higher except for the Yumizen G PT Liq 4 lots 920611 & 920705. For the other lots, including 920611 & 920705, please refer to the previous revision of values sheet.



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 1 0 4 3 4 5 6 4 1 2 0 4 2 3 D



D 5 1 4 7 2 2 1 0 3 5 0 5 9 1 2 9 D



D 6 1 9 3 1 4 3 2 1 4 0 3 2 0 5 3 D



D 7 1 2 9 1 9 4 1 0 8 1 6 2 0 4 1 D



D 8 0 6 8 1 0 6 1 5 9 4 4 8 6 0 7 D



D 9 4 3 2 5 8 4 0 8 0 1 9 0 0 0 0 D



D 0 9 9 9 3 6 4 D

Reference No.:		1300036412	
Lot No:		910434 (rev. 3)	
Expiry date:		2023-04-30	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	14,7-22,1
	Prothrombin %	%	35-59
	Prothrombin INR	INR	1,29-1,93
PT Liq (*)	Prothrombin Time	sec	14,3-21,4 (*)
	Prothrombin %	%	32-53 (*)
	Prothrombin INR	INR	1,29-1,94 (*)
PT R	Prothrombin Time	sec	10,8-16,2
	Prothrombin %	%	41-68
	Prothrombin INR	INR	1,06-1,59
APTT	Partial Thromboplastin Time	sec	44,8-60,7
APTT Liq	Partial Thromboplastin Time	sec	43,2-58,4
FIB Cl	Fibrinogen	g/L	0,8-1,9
TT	Thrombin Time	sec	NA
AT	Antithrombin	%	44-73

Important notice

The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.

The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.

(*) Warning! This range has to be used only with Yumizen G PT Liq 4 lots : 920315 and higher except for the Yumizen G PT Liq 4 lots 920611 & 920705. For the other lots, including 920611 & 920705, please refer to the previous revision of values sheet.