



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4 9 1 1 2 2 5 4 6 4 1 2 1 2 2 3 D



D 5 1 1 9 1 6 1 0 6 2 1 0 4 0 8 7 D



D 6 1 3 1 1 2 4 1 6 8 0 6 3 1 0 6 D



D 7 0 8 8 1 3 2 0 9 2 1 2 4 0 6 9 D



D 8 1 1 5 0 8 3 1 2 4 3 3 6 4 5 5 D



D 9 3 2 2 4 3 6 2 0 0 3 0 0 1 8 6 D



D 0 2 7 8 3 3 7 D

Reference No.:		1300036412	
Lot No:		911225 (rev. 4)	
Expiry date:		2023-12-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	11,9-16,1
	Prothrombin %	%	62-104
	Prothrombin INR	INR	0,87-1,31
PT Liq	Prothrombin Time	sec	12,4-16,8
	Prothrombin %	%	63-106
	Prothrombin INR	INR	0,88-1,32
PT R	Prothrombin Time	sec	9,2-12,4
	Prothrombin %	%	69-115
	Prothrombin INR	INR	0,83-1,24
APTT	Partial Thromboplastin Time	sec	33,6-45,5
APTT Liq	Partial Thromboplastin Time	sec	32,2-43,6
FIB Cl	Fibrinogen	g/L	2,0-3,0
TT	Thrombin Time	sec	18,6-27,8
AT	Antithrombin	%	86-116
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 1 1 2 2 5 5 6 4 1 2 1 2 2 3 D



D 5 1 4 3 2 1 4 0 4 1 0 6 9 1 1 4 D



D 6 1 7 1 1 6 7 2 5 1 0 3 3 0 5 4 D



D 7 1 4 6 2 1 9 1 1 1 1 6 7 0 4 2 D



D 8 0 7 0 1 0 6 1 5 9 5 1 7 7 0 0 D



D 9 4 9 6 6 7 1 0 8 0 1 9 0 0 0 0 D



D 0 9 9 9 3 4 5 D

Reference No.:		1300036412	
Lot No:		911225 (rev. 4)	
Expiry date:		2023-12-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	14,3-21,4
	Prothrombin %	%	41-69
	Prothrombin INR	INR	1,14-1,71
PT Liq	Prothrombin Time	sec	16,7-25,1
	Prothrombin %	%	33-54
	Prothrombin INR	INR	1,46-2,19
PT R	Prothrombin Time	sec	11,1-16,7
	Prothrombin %	%	42-70
	Prothrombin INR	INR	1,06-1,59
APTT	Partial Thromboplastin Time	sec	51,7-70,0
APTT Liq	Partial Thromboplastin Time	sec	49,6-67,1
FIB Cl	Fibrinogen	g/L	0,8-1,9
TT	Thrombin Time	sec	NA
AT	Antithrombin	%	37-62
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			