

ABX Minocal



LOT

CX 472

CAL



(Exp.) 2022-08-05

(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITES UNITS | ABX Lyse | | | | | | | | TOLERANCES TOLERANCE | |
|--------------------------|---|-----------------|--------------|--------|--------|-----------------------------|--------|--------|--------|-------------------------|---------|
| | | MICROS | MICROS | MICROS | MICROS | MICROS | PENTRA | PENTRA | PENTRA | | PENTRA |
| | | ES60 Care ST | CRP 200 | | | 60 8P 60 18P ADVIA 60 | | | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 9.2 | 8.8 | | | 9.0 | | | | | ± 0.2 |
| GR RBC | 10 ⁹ /mm ³ ;10 ¹² /l | 4.51 | 4.50 | | | 4.50 | | | | | ± 0.06 |
| HB HGB | g/dl | 13.6 | 13.7 | | | 13.5 | | | | | ± 0.2 |
| | g/l | 136 | 137 | | | 135 | | | | | ± 2 |
| | mmol/l | 8.45 | 8.51 | | | 8.38 | | | | | ± 0.12 |
| HT HCT | % | 38.6 | 38.9 | | | 38.7 | | | | | ± 1.0 |
| | l/l | 0.386 | 0.389 | | | 0.387 | | | | | ± 0.010 |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 251 | 248 | | | 245 | | | | | ± 10 |
| VMP MPV | μm ³ ; fl | 9.0 | 8.5 | | | 7.9 | | | | | ± 0.5 |

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308