

ABX Minocal

LOT		CX 473		CAL		(Exp.) 2022-09-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.1			9.0					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.45			4.45					± 0.06		
HB HGB	g/dl	13.2	13.4			13.3					± 0.2		
	g/l	132	134			133					± 2		
	mmol/l	8.20	8.32			8.26					± 0.12		
HT HCT	%	37.5	38.0			37.8					± 1.0		
	l/l	0.375	0.380			0.378					± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	251			250					± 10		
VMP MPV	µm ³ ; fl	8.2	8.2			7.4					± 0.5		

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308