

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)					
Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.0			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.55			4.51					± 0.06
HB HGB	g/dl	13.4	13.3			13.3					± 0.2
	g/l	134	133			133					± 2
	mmol/l	8.32	8.26			8.26					± 0.12
HT HCT	%	38.2	37.8			37.7					± 1.0
	l/l	0.382	0.378			0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240	245			237					± 10
VMP MPV	µm ³ ; fl	9.0	8.9			8.3					± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48